

## 1. PRODUCT IDENTIFICATION

Premium Kohno TG dental implants are implantable devices for the rehabilitation of patients affected by total or partial edentulism. The implants are intended to be inserted into the mandibular or maxillary bone (that is, as implant fixtures) via a surgical operation.

The fixtures have a connection in the crown part to hold the abutment implant which acts as a support for the dental prosthesis. The purpose of dental prostheses is to restore the dental function to patients from an aesthetic, phonetic and masticatory point of view.

In implantoprosthesis rehabilitation using Kohno TG implants, exclusively original Sweden & Martina prosthetic components must be used. Use of components that are not produced by Sweden & Martina limits their liability and renders the product Warranty void (see the later section "Liability for Defective Products and Terms of Warranty").

For surgical insertion of the fixture, appropriate surgical instruments must be used, available either singularly or in kits. Use of original surgical accessories manufactured by Sweden & Martina is recommended. Sweden & Martina does not have any liability if non-original instruments are used.

Kohno TG implants can be inserted in different sites in the oral cavity employing various techniques and then connected to the prosthesis at a later time. Implants (or, more correctly, implant bodies or fixtures) a truncated conical form. They are screw-shaped with a transgingival coronal part and a connection for the prosthetic superstructure featuring a Morse taper with a hexagon below. Depending on operating time (for functionality restoration), they can be rehabilitated with immediate, anticipated or deferred loading. Kohno TG implants can be inserted in sites that are already edentulous or in post-extraction sites, both immediately (insertion of the implant takes place at the same time as tooth or root is removed) or at a deferred time (a period of about 3 weeks is normally allowed to pass between extraction and insertion of the implant fixture).

## 2. INTENDED USE

Kohno TG implant fixtures are medical devices intended for long term implantations. All the fixtures are sold in single-use, sterile packaging. The purpose of the fixtures is to replace the missing roots of the tooth.

All the fixtures are sold complete with their respective cover screws (also called surgical screws). The cover screws are also implantable medical devices of surgical type, designed to remain in the oral cavity for a duration that can exceed 30 days. Cover screws are also available in individual packs. In this case also the pack is sterile.

Sweden & Martina declares that it is the manufacturer of the Kohno TG implants and attributes the risk classes given in Table 01. The dental implants, although intended for implantation in all subjects who satisfy the appropriate therapeutic indications, must be used exclusively by professional, medically-qualified personnel having the necessary qualifications and approvals.

## 3. MANUFACTURER'S DETAILS

The Manufacturer of the Kohno TG implant fixtures is:

**Sweden & Martina S.p.A.**  
Via Veneto 10 - 35020 Due Carrare (Padova) - Italia  
Tel. +39 049.91.24.300 - Fax +39 049.91.24.290  
e-mail: info@sweden-martina.com - www.sweden-martina.com

## 4. RAW MATERIALS USED

The materials used for manufacturing the Kohno TG dental implants were selected on the basis of the properties required for their intended use, in accordance with Regulation (EU) 2017/745.

The implants are produced in grade 4 titanium, according to the harmonised normatives.

Allergies to titanium are very rare but nevertheless possible. For this reason, it is always necessary to check in advance with patients that they do not suffer from this type of allergy.

## 5. DESCRIPTION

The information contained in these Instructions for Use complete the indications given in the catalogues and manuals. If you do not have copies, they can be requested from Sweden & Martina S.p.A.

Kohno TG implants are screw implants with a number of features specially studied to simplify the surgical procedure as in the most recent implantology protocols.

The way the implants are designed (see Fig. 01) guarantees a great primary stability, even in the presence of poor quality bone, thanks to the presence of a marked spiral in both the body and in the apical portion. The apical incisions guarantee an easier insertion into dense bone, making the insertion procedure extremely simple and non-traumatic. Implants are available, on the basis of single specifications in diameters of 3.30, 3.80, 4.25 and 5.00 mm and in heights varying from 6 to 15 mm. The complete range can be seen in the related catalogue.

The thread has a pitch of 1 mm, a tapering profile in the apical direction, convex in the crown direction, and joined to the implant body by a circular section. This profile avoids the bone to be subjected to trauma when the load is applied, creating the perfect conditions for complete bone integration.

The apical incisions allow the bone to be tapped, providing three decompression and release areas for the bone fragments, and improve primary stability, helping also prevent rotation of the implant during the second surgical phase while performing the manoeuvres for screwing and unscrewing the components connected to it.

Pre-tapping of the bone is always appropriate in the case of very dense bone (D1).

While the two types of implantation morphologies offer different surgical options as a result of the differences in their morphologies, they nevertheless share a common prosthetic platform, and this therefore makes the work of the prosthodontist easier.

The implants have a conical internal connection, anchored by means of an internal hexagon of ample dimensions at the base of the cone, that permits the posts to be repositioned at the angles required.

The implants of diameters 3.80 and 4.25 all have the same connection platform, and adopt the same prosthetic solutions. The diameter of the platform is 4.80 mm, while the hexagon underneath the Morse taper measures 2.3 mm. The fixing screw is M 1.8.

Implants of 5.00 mm diameter have a 5.80 mm platform and therefore need dedicated posts different from those used for implants of smaller diameters, that in this case would offer too fine a profile with respect to the total dimensions of the platform. The hexagon of these implants is 2.5 mm in size. The fixing screw is M 2.

The special, patented, design of the connection allows the head of the screw used to fix the post to descend inside the internal hole of the platform, minimising in this way the stress it is subjected to, with significant benefits:

- This feature combined with the presence of the Morse taper means that the risks of unscrewing and breakage are minimal while, at the same time, the head of the fixing screw is no longer an obstacle for the prosthodontist during the preparation of the abutment.

Within this scheme, fixtures of reduced height (6, 7 and 8.5 mm) are available to complete the implantoprosthesis solution, and can be used in accordance with the most recent clinical protocols in all cases where the vertical bone is of reduced dimension. The slight tapering of the apex makes insertion of the fixture easier, while the thread pitch and depth guarantee optimum primary stability. Given the reduced dimension of these implants, they should only be used for supporting multiple prostheses, together with implants of greater dimensions.

The measurement of the implant diameter always refers to the spiral diameter measured at the widest point of the spiral. The size of a prosthetic platform is calculated as the diameter of the emergence at its widest point, at the perimeter of the bevel constituting the closing shoulder of the prosthetic posts. By implant length is always intended the length of the fixture referring exclusively to the endosseous part. The total length is given by the endosseous length to which is added the 2.2 mm of the emerging part.

Kohno TG implants are available with Zirti (sandblasted and acid etched body) surface treatment and UTM (Ultrathin Threaded Microsurface) neck.

## 6. METHOD OF USE

Modern implantology, whether using immediate or deferred loading, is largely an experimented reliable discipline able to resolve almost all problems related to edentulism, whether they be functional or aesthetic in nature.

Kohno TG implants have been specially designed for one-step surgical procedures: the implant is inserted into the bone with the coronal part emerging. Thus the implant can be left to heal for the normal healing period and bone integration, or loaded immediately with an appropriate dental post, provisional or definitive, as needed.

Implants are inserted into the bone following surgical protocols that must take into account the quantity and quality of the receiving bone, the implant type, and the possible need for regenerative therapies. A site is created in the patient's bone (corresponding to the site for the tooth to be replaced or built anew altogether) using a series of calibrated bone millers or appropriate instruments such as bone-expanders, bone compactors etc. In order for the implant to osseointegrate, a good primary stability is required with little or no movement - if movement is present it must not exceed a few microns. The bone-implant interface is therefore to the order of mill-microns, otherwise the implant risks being fibrointegrated with no osseointegration.

In general, the masticatory load with the fixed prosthesis takes place in a second phase, after 2-3 months for the mandibular bone and 4-6 months for the upper maxillary bone. In certain cases, immediate loading of the implant is possible, however certain fundamental criteria must be met:

- presence of a certain amount of bone
- primary stability of implants after insertion
- good parodontal (gingival) support
- absence of bruxism (tooth greying) or serious malocclusion
- presence of a good occlusal balance (correct masticatory occlusal plane).

A serious assessment is therefore required of the specialist who, carrying out the necessary examinations with the help of appropriate instruments, must verify the coexistence of all these factors. If not, the choice must fall on "traditional" techniques (of "submerged" or "non-submerged" type), that is, using implants that require a longer waiting time but that are stronger for supporting the masticatory load.

Implants can replace a single tooth (crown on implant), a group of teeth close together (bridge on implant) or a whole dental arch, or they can be used to stabilise a full upper or lower overdenture prosthesis.

The clinical diagnosis for determining the type of implant and its size depend on: the site for which the implant is intended; the anatomy of the receiving bone; the number of implants; the technically-motivated choice of the most appropriate protocol from those indicated above. The choice must be made exclusively by the doctor performing the operation, who must have suitable

preparation, and programme in advance the appropriate prosthetic rehabilitations. Whenever possible, implants of the largest diameter for the crest thickness must be used. Bear in mind the usage limitations for short implants (heights 6, 7 and 8.5 mm) and narrow implants (diameter 3.30 mm), indicated in the previous paragraph.

## 6.1 Pre-surgery Planning and Preparation

During the phase preceding the intervention, the following is required:

- General medical and dental history, general medical examination, clinical examinations (full blood tests) and radiological examinations, TAC and consultation with family doctor
- Patient information (indications, counter-indications, clinical situation, expectations, normal percentages of success and failure, necessity for periodical post-surgery check-ups)
- Hygiene programme, with periodical interventions (if any)
- Adoption of the necessary pharmacological prescriptions
- Pre-prosthetic surgical planning in collaboration with dental technician
- Assessment of risks relating to inadequate treatments of soft and hard tissues
- Choice of anaesthetic and sedation techniques, and amount of monitoring necessary
- Prosthetic planning in collaboration with dental technician.

## 6.2 Surgical Intervention

The surgical techniques for implants are taught in University establishments to students who are graduating in dentistry. Nevertheless, the following factors must be born in mind:

- Hard and soft tissue must be treated with extreme care, taking all the necessary precautions in order to obtain a good integration of the implant.
- The normal biological principles for osseointegration must be respected.
- Thermal traumas must be avoided, they could cause necrosis and reduce the possibility of osseointegration. For this reason, a reduced drilling speed must be appropriate and milling tools with cutters in optimum conditions. Drilling must be carried out intermittently, cooling the site adequately with the irrigation necessary, and widening the hole using milling tools with gradually increasing cooling diameters.
- It is advisable to gather and keep on file a complete clinical, radiological and x-ray documentation
- It is indispensable to respect the healing times recommended for implant surgery and to check periodically - also by means of radiographic controls - the state of progression of osseointegration.

## 6.3 Instructions relating to Product Handling and Storage

Implantological interventions must be carried out in a suitable, adequately aseptic environment. It is recommended to always cover surfaces with sterile drapes, to cover the dental unit and the micromotor with appropriate coverings to isolate the operating field covering the patient with appropriate surgical gowns, to wear sterile gloves and open the sterile packages containing the instruments only just before their use.

Kohno TG implants are packaged in sterile vials, placed inside a blister with a Tyvek seal. This blister is in turn contained in a little box that forms the external enclosure. The package also contains the adhesive labels for use on the patient's medical card. The blister preserves the sterile conditions and is shaped and pre-formed in a way to prevent the vial from moving as much as possible, while at the same time making it easy to extract. The blister is sealed with a Tyvek film (Figs. 02 and 03). The sterile blister must only be opened in controlled asepsis conditions. Remove the vial from its seat. The vials containing the implants must be opened only in a sterile environment, immediately before insertion of the fixtures into their sites. Inside the sterile vials, the fixtures are enclosed in special sleeves in titanium keeping them straight, with the connection visible, ready to be taken up by the surgical instruments.

Kohno TG implants are designed for use in a mouthless surgical procedure: special drivers (Sweden & Martina SpA Easy Insert) engage directly the internal hexagons of the connection, allowing the fixtures to be extracted from the phials without the need for contact by hand or other instruments, so avoiding any risk of contamination before use. The drivers have been specially studied to avoid problems relating to deformation of the connections, or excess engagement during the surgical phase, so as to limit mechanical damage.

Any contact, even accidental, with the surfaces of the implant before its insertion into the surgical site nullifies the ideal surface conditions produced by the surface treatment procedure.

Should it be necessary to manipulate the implant to insert it into the prepared site, you must use only clean and sterilised tweezers in titanium. All contacts between the surfaces of the implant and the epithelial and connective tissue must be avoided, because they could prejudice the success of the intervention.

At the end of the intervention, if the implant is to be submerged, before closing the flaps, the connection must be sealed well using the proper cap screw. The cap screw is located in a special location, inside the light blue cap (Fig. 04) used to close the vial. A small label indicates its presence. The cap screw can be extracted using an appropriate electric screw driver and carried directly into the implant. At the end of the intervention, the flaps must be repositioned and closed. Suturing should be performed as normal.

Every package displays the code and a contents description, the batch number, the "sterile" label and its expiry date. This same data is also indicated on the labels for use on the patient's medical card and must always be stated by the doctor in any communication in its relation.

The packaging is in compliance with European regulations.

Implants must be stored in a cool dry place, protected from direct sunlight, water and heat sources.

## 7. STERILISATION

Kohno TG implants are sterilised using Beta rays. The expiry date is given on the packaging. The sterile blister must be opened only at the moment it has to be used in the operation. Before opening, make sure that the packaging is perfectly intact. Any damage could compromise the sterility of the implant, and therefore the success of the intervention. Implants that have already been used, or that are not in a sterile condition, must never be reused.

The device is for single-use only: its reuse is not permitted as it could lead to loss of the implant and cross-infections.

On the bottom of the vial there is a round label (or sticker). This label certifies that sterilisation has been obtained through radiation. In fact, this label is yellow in the beginning and turns red under the effect of the radiation, thereby confirming that sterilisation has taken place.

## 8. COUNTER-INDICATIONS

Insertion of implants and prosthesis implants is counter-indicated in patients presenting a poor general health condition, scarce or inadequate oral hygiene, or where it is impossible or difficult to monitor their general conditions, or in patients who have previously been subjected to organ transplants. Patients with psychiatric problems must also be excluded, as well as those prone to alcohol or drugs abuse, who are little motivated or not sufficiently co-operative. Patients whose gums are in a bad condition must be treated and their condition recuperated in advance. In cases where the receiving bone contains insufficient material or is of such a poor quality that the implant stability could be jeopardised, an appropriate guided regeneration of the tissue must be carried out in advance. Other counter-indications include: allergies to titanium, acute or chronic infective diseases, chronic sub-acute maxillary osteitis, systemic diseases, endocrine disorders, diseases leading to microvascular disorders, pregnancy, breastfeeding, previous exposures to radiation, haemophilia, granulocytopenia, use of steroids, diabetes mellitus, renal insufficiency, bone fibrous dysplasia. All the normal counter-indications for oral surgical operations must also be taken into account. Patients must not be subjected to interventions if they are undergoing anticoagulant, anticonvulsant or immunosuppressive therapies, if inflammatory-infectious processes are present in the oral cavity, or if their creatinine or BUN values are outside the normal range. Other patients who must be excluded are those with cardiovascular disease, hypertension, thyroid or parathyroid diseases, malignant tumours encountered within the 5 years preceding intervention, or node enlargements.

Chemotherapies reduce or counter osseointegration capability, therefore patients undergoing such treatments must be accurately evaluated before implantoprosthesis rehabilitations. In the literature, numerous cases of periimplant osteonecrosis were reported in patients that had been administered with bisphosphonates, in particular in the lower jaw and mainly in case of intravenous administration.

Failures of implants inserted in sites where a root canal therapy have been previously performed, are reported in the literature. Planning implant placements, any previous endodontic therapies should be carefully evaluated in the patient's medical history. Unexpected implant failures have been reported in the literature related to patients who regularly assume, or even only for repeated periods, proton pump inhibitors. Planning implant-prosthetic rehabilitations, carefully evaluate the possible assumption of such drugs by the patient.

## 9. SPECIAL WARNINGS

As a precaution, after the intervention, the patient must avoid activities requiring physical effort.

When tightening the cap screws, post screws or prosthetic screws, you must adhere strictly to the tightening torque recommended in the related Instructions for Use. Too high a tightening torque could weaken the mechanical structure of the screw and compromise the prosthetic stability, causing possible damage to the implant connection.

## 10. SECONDARY EFFECTS

After dental implant operations the following could occur: bone crest loss, permanent numbness, dysaesthesia, local or systematic infections, exfoliation, hyperplasia, oroantral and oronasal fistula. Temporary complications can occur such as pain, swelling, pronunciation problems and gingivitis. The risks related to an implantological intervention include: perforation of the lip or tongue plate, bone fractures, implant fractures, fractures of the upper structures, aesthetic problems, accidental perforation of sinus, nerve damage and compromise of natural dentition. The following pathophysiological complications can increase the degree of risk: cardiovascular insufficiency, coronary disorders, arrhythmia, chronic respiratory or lung diseases, gastrointestinal diseases, hepatitis, intestinal inflammations, chronic kidney insufficiency and urinary system disorders, endocrine disorders, diabetes, thyroid diseases, blood disorders, anaemia, leukaemia, coagulation disorders, osteoporosis or musculoskeletal arthritis, heart attacks, neurological disorders, mental retardation, paralysis.

## 11. MAINTENANCE

Complications linked to implant prostheses are documented in the related literature. These complications can lead to a lack of osseointegration and a failure of the implant. A correct up-keep on the part of the patient, with a regular attention to dental hygiene at home, combined with periodic check-ups and visits to a professional hygienist lengthen the useful life of the device. Complications such as, for example, loosening of the screws securing the prosthesis to the implant, or bone re-absorption causing loss of mucosa support for removable prostheses, can easily be detected by regular control visits.

Should it be necessary to tighten the abutment or the prosthetic screws, these operations must be carried out by the doctor using the appropriate devices that are able to verify the tightening torque. Devices must be calibrated on a regular basis. Should the patient become aware that any of the conditions above have occurred, they should contact their doctor as soon as possible so that the prosthesis can be restored to its proper functional condition. Any delay in requesting medical intervention could, in the first instance, lead to the fastening screw or the prosthesis fracturing and, in the second instance, to loss of the implant, affecting the rehabilitation result achieved. Doctors must therefore prepare patients for these circumstances. Complications can be biological (integration loss) or mechanical (component fracture due to excessive load). If no complications occur, the duration of the devices and of the prosthesis as a whole depends on the mechanical resistance of the device with respect to the accumulated fatigue.

Sweden & Martina has subjected Kohno TG implants to the required fatigue resistance tests at 5,000,000 cycles, and the implants passed the test.

Fixture tests are carried out in compliance with the specific normative and the results assessed by performing calculations on the finished elements.

#### 12. EXPIRY DATE

Implants must not be used after the expiry date indicated.

#### 13. LEGAL REFERENCES

The design and production of Kohno TG implant fixtures are carried out in conformity with the directives and most up-to-date harmonised normatives regarding the materials used, production processes, sterilisation, information provided and packaging.

#### 14. WASTE DISPOSAL PROCEDURES

Fixture implants, if removed from the oral cavity as a result of a biological or mechanical failure, must be treated as organic waste for their disposal, according to the laws that apply locally.

On the other hand, if the implants are sent to Sweden & Martina with a request for execution of a Surf Test, the protocol given on the website [www.sweden-martina.com](http://www.sweden-martina.com) must be followed.

#### 15. LIABILITY FOR DEFECTIVE PRODUCTS AND TERMS OF WARRANTY

Excellent patient care and attention to their needs are necessary conditions for the success of the implant. It is therefore necessary to select the patient carefully, inform them of the inherent risks and of the duties associated with the treatment, encouraging them to co-operate with the dentist to achieve a good outcome from the treatment. The patient must therefore maintain a good level of oral hygiene - confirmed by means of regular check-ups and control visits - that must be guaranteed and documented, along with the pre- and post-surgical directions and prescriptions.

The instructions provided by Sweden & Martina are available at the moment of treatment and have been accepted by the Dental Practice. These instructions must be observed and applied during all the care phases: from the patient medical history stage to the post-surgery check-ups.

The Warranty covers exclusively defects that are established to be attributed production-related and on submission of the piece identified by item and batch code, within validity period of the Warranty. The Warranty Conditions are available on the [www.sweden-martina.com](http://www.sweden-martina.com) website.

#### 16. DATE AND VALIDITY OF INSTRUCTIONS FOR USE

These Instructions for Use have validity and effect from the month of July 2021.

picture 01



picture 02



picture 03



picture 04




Table 01

Device	Packaging	Regulation (EU) 2017/745	Classification Rule	Risk Class
Implant fixtures for dental use, belonging to the Premium-Kohno TG implantological system	Sterile, single-use packages. Fixtures come complete with cover screw	Implantable devices intended for long term use (over 30 days)	8	Iib
Cover screws	Sold either complete with appropriate fittings or individually (single-use, sterile package)	Implantable devices intended for long term use (over 30 days)	8	Iib

Table 02

ø 3.80 TG	ø 4.25 TG	ø 5.00 TG
Green	Blue	Magenta

EXPLANATION OF SYMBOLS		
	Caution! See instruction for use	✓
	Batch number	✓
	Code	✓
	Manufacturer	✓
	Country of manufacture	✓
	UDI code, Unique Device Identification	✓
	Medical Device	✓
	Consult instruction for use <a href="http://www.sweden-martina.com">www.sweden-martina.com</a>	✓
	CE marking Where applicable: The identification number of the Notified Body shall follow this symbol.	✓
	American federal law restricts this device to sale by or by order of a professional practitioner	✓
	Do not resterilize	✓
	Disposable product, do not reuse	✓
	Do not use if the packaging is damaged	✓
	Sterilized with ionizing radiation	✓
	Single sterile barrier system with protective packaging inside	✓
	Expiry date after which the product must not be used	✓



Sweden & Martina S.p.A.  
Via Veneto, 10 - 35020 Due Carrare (Padova) - Italy  
Tel. +39.049.9124300 - Fax + 39.049.9124290  
e-mail: [info@sweden-martina.com](mailto:info@sweden-martina.com)  
[www.sweden-martina.com](http://www.sweden-martina.com)