

Surgical manual

PRAMA



Prama



The implants	4
Clinical indications for resorting to implantoprosthesis therapies	4
Side and secondary effects	5
General indications	6
Method of use	7
Prama implants	8
Prama RF implants	9
Surface	10
Key to the implant codes	11
Colour codes	11



The range	12
Prama	12
Prama RF	13



Surgical instruments	14
Surgical kit	14
General indications	18
Drills	19
Precision drill FS-230	20
Pilot drill FPT3-200-LXS	20
Intermediate drills	22
Other intermediate drills	22
Final cylindrical drills and related stops	24
Final conical drills and related stops	26
Reply: replies for Prama RF implants	30
Countersink drills	32
Drills for distal sectors	34
Bone taps	35
Osteotomes	36
Parallelism pins	38
Depth gauge PROF3	39
X-ray templates	39
Easy Insert Driver	40
Drivers for fixation screws	42
Dynamometric ratchet CRI5-KIT	46
Torque wrench with control lever TWL	48
Cleaning, disinfection, sterilisation and storage of the kit and of the surgical instruments	49
Shorty Drilling Kit	54
Drilling Kit Syra Short	58



Surgical sequences	62
Preparation of the implant site	62
Surgical sequences	62
Surgical sequence for Prama implants	64
Surgical sequence for Prama RF implants	76



Surgical procedures	88
Implant insertion	88
Intraoperative removal of the implants if necessary	90



General indications	91
Maintenance of the prosthesis	91
Responsibility for defective products and warranty terms	91
Disposal	91
Material composition	92
Identification of the manufacturer	93



Bibliography on Sweden & Martina implants since 2013	96
--	----

Clinical indications for resorting to implantoprosthesis

When assessing the patient, in addition to his/her eligibility as regards implant-prosthetic rehabilitation, it is usually necessary to consider the contraindications that apply to oral surgery procedures in general. These include:

- clotting disorders, anticoagulant therapy;
- healing or bone regeneration disorders;
- decompensated diabetes mellitus;
- metabolic or systemic diseases that compromise tissue regeneration with a particular influence on healing and bone regeneration;
- alcohol abuse, smoking and use of drugs;
- immunosuppressive therapy, such as: chemotherapy and radiotherapy;
- infections and inflammations, such as periodontitis and gingivitis;
- poor oral hygiene;
- inadequate motivation;
- occlusion and/or articulation disorders as well as an inadequate interocclusal space;
- inadequate alveolar process

It is contraindicated to fit implants and implant restorations in patients with poor general or oral health, those who are unable to monitor their general conditions properly or those who have had organ transplants. Psychologically unstable patients, alcohol or drug abusers, and poorly motivated or uncooperative patients should also be considered unsuitable for this kind of treatment. Patients with poor periodontal health should first be treated and allowed to recover. In the presence of a lack of bone substance or poor quality of the receiving bone, such as to compromise the stability of the implant, suitable guided tissue regeneration must be performed prior to implant treatment.

Contraindications also include: bruxism, allergy to titanium (extremely rare), acute or chronic infectious diseases, sub-acute chronic maxillary osteitis, systemic diseases, endocrine disorders, diseases resulting in microvascular disorders, pregnancy, breastfeeding, previous exposure to radiation, haemophilia, neutropenia, steroid use, diabetes mellitus, kidney failure and fibrous dysplasia.

The normal contraindications common to all oral surgery must also be observed. Surgery is not recommended for patients on anti-coagulant, anticonvulsant and immunosuppressant therapies, with active inflammatory-infective processes of the oral cavity, and patients with BUN and creatinine values outside the norm. Patients with cardiovascular disease, hypertension, thyroid or parathyroid diseases, malignant tumours found in the 5 years preceding the operation, or nodular swellings must also be rejected. Chemotherapies reduce or eliminate the ability of osseointegration, therefore patients undergoing these treatments must be carefully screened before being rehabilitated with oral implantoprosthesis. Numerous cases of bisphosphonate-associated periimplant osteonecrosis of the mandible have been reported in the literature. This problem particularly applies to patients treated intravenously. As a post-operative precaution, the patient must avoid any kind of strenuous physical activity.

Side and secondary effects

Situations that may occur after surgical procedures include temporary local swelling, oedema, haematoma, temporary sensitivity alterations, temporary masticatory limitations, post-surgical micro-hemorrhages in the following 12-24 hours. The patient may also experience pain, speech problems, gingivitis, loss of bone crest, permanent paresthesia, dysesthesia, local or systemic infections, exfoliation, hyperplasia, and oronasal and oroantral fistulas, perforation of the labial or lingual plate, perforation of the Schneider membrane, bone fractures, implant fractures, fractures of the over-structures, aesthetic problems, unnoticed perforation of the nasal sinus, nerve injuries, impairment of natural dentition. The following pathophysiological problems can increase the risks: cardiovascular failure, coronary disease, arrhythmia, pulmonary or chronic respiratory disease, gastrointestinal disease, hepatitis, inflammatory bowel disease, chronic kidney failure and disorders of the urinary system, endocrine disorders, diabetes, thyroid diseases, hematologic disorders, anaemia, leukaemia, coagulation problems, osteoporosis or musculoskeletal arthritis, stroke, neurological disorders, mental retardation, paralysis.


Before proceeding, it is important to perform a careful pre-operative analysis of the patient's medical history to verify his or her suitability for the implant treatment. It is also recommended to collect and file all the clinical, radiological and radiographic records. After making models of the two arches, the best position and orientation of the chosen implants will be evaluated based on the occlusal plane and on a correct distribution of the forces. In this phase, a surgical stent may be created to guide the specialist to correctly position the implants during the operation. Depending on the specific case, a decision will be made on whether to use a single or double phase surgical procedure, using titanium cylinders (code DIM) to make the radiological/surgical stent.




A radiological and surgical stent can be made by using the special cylinders in titanium (code DIM), which can be used to obtain an ideal positioning of the implants in terms of biomechanics and aesthetics.

In addition to an oral examination, both clinical and with x-rays, it is recommended to take a T.C. scan of the interested area; once the x-rays and scans have been obtained, the specialist can identify the most suitable implant with the help of convenient transparent radiographic guides.


The pre-operative study of the T.C. Dentscan allows identifying the type of bone present in the insertion point of the implant. The choice of the surgical procedure must take into consideration the type of bone present. The bone is normally classified into 4 types according to the density. The classification (according to Carl Misch) is the following:




BONE D1: all cortical bone.



BONE D2: a core of bone marrow enclosed in a shell of cortical bone.



BONE D3: all bone marrow without crest cortical.




BONE D4: all bone marrow with very poor mineralisation.

General indications

Prama implant fixtures mesh are long-term implantable medical devices. All the fixtures are sold in single-use sterile packs. The function of the fixtures is to replace missing dental roots. The fixtures have a connection in the crown part for receiving an implant post aimed at supporting a dental prosthesis. In implant-prosthetic rehabilitation with Prama implants, exclusively original prosthetic components by Sweden & Martina must be used.

Use of non-original components limits the responsibility of Sweden & Martina S.p.A. and renders the product warranty void. The implants have a cylindrical shape (Prama) and conical shape (Prama RF), they are screw shaped with an external thread and have a hexagonal internal connection for connecting the prosthetic components. Prama implants can be inserted in both edentulous and post-extraction sites, either immediate (insertion of the implant at the same time as the removal of the tooth or root), or deferred (normally about 3 weeks between extraction and insertion of the implant fixture).

All the fixtures are sold with the respective closing cover screws (also called, surgical cover screws). The surgical cover screws are also medical devices that can be implanted surgically. They are designed to remain in the oral cavity for more than 30 days. The surgical cover screws can also be sold individually.



Method of use

The method to be used is the one stage, surgical technique: insertion of the implant, closure of the connection with a healing abutments, surgical cover screw. Alternatively, in the presence of suitable therapeutic indications, it can be loaded immediately with an appropriate temporary or permanent dental post, depending on the case.

Implants are inserted in the bone based on surgical protocols that must be considered according to the quantity and quality of the receiving bone, the implant, and the possible need for regenerative therapies. The “implantologist” or dental surgeon creates a site in the patient’s bone (corresponding to the new tooth to be placed or replaced), by using a series of calibrated drills or suitable instruments such as bone expanders, bone compactors or similar instruments.

The necessary conditions for the success of the implant are:

- the presence of a certain amount of bone;
- good periodontal (gingival) support;
- no bruxism (teeth grinding) or serious malocclusion;
- the presence of good occlusal balance (correct masticatory occlusal plane).

Prima implants have been tested in a wide range of clinical situations:

- standard operating procedures involving the double or single surgical phase;
- immediate and early loading;
- post-extraction situations, even combined with immediate loading.

Generally, masticatory loading with a fixed prosthesis occurs at a second stage, after 2 to 3 months for the mandible and after 4 to 6 months for the upper jaw. In some cases, but not all, immediate loading of the implants is possible; to do this it requires good primary stability, with no mobility or movement limited to a few microns. The bone-implant interface must therefore be of the order of a few millimicrons, otherwise there is the risk of fibrous integration.



In the case of temporary single crowns with immediate loading, it is recommended to avoid a direct occlusion with the antagonist; in the case of multiple solutions the temporary prosthesis has to be splinted in a single structure.

The clinical indication for choosing the Prima implant depends on the site in which the implant is to be inserted, on the anatomy of the receiving bone and on the technique chosen from among those mentioned above. The choice must be made exclusively by the doctor, who must have the suitable training and experience and must plan the prosthetic rehabilitations beforehand.


Sweden & Martina has conducted 5.000.000-cycle fatigue resistance tests on Prima implants.

The implants passed the test. Fatigue tests are conducted according to the standards and evaluated further with finite element calculations.


Prama implants



The convergent neck of the Prama implants is characterized by a first cylindrical section 0.80 mm high, followed by a section with hyperbolic geometry. The connection diameter is of 3.40 mm in all the available implant diameters.



The Prama implant is characterized by the Collex connection, with internal hexagon and collar for the prosthetic support which gives resilience and stability to the prosthesis and acts as a guide and connection for the Easy Insert drivers.




The spire of the Prama implant has an asymmetric profile and its thread has a pitch of 1.00 mm and a depth of 0.40 mm.




The apex of Prama implants has three incisions that increase its penetration capacity and guarantee a good self-threading capacity, improve the primary stability and offers three zones to decompress the clot.

Prama RF implants




The neck of the Prama RF implants has the same characteristics of the Prama one: the convergent section has with hyperbolic geometry and is followed by a cylindrical section 0.80 mm high.



The Prama RF implants, as well as the Prama ones, are characterized by the Collex connection, with internal hexagon and a collar for the prosthetic support.



The thread of the Prama RF implants is characterized by a triangular profile, a pitch of 1.50 mm and a depth of 0.40 mm.





The apex of Prama RF implants has two incisions that increase its penetration capacity and non-rotational property. The hemispherical apex makes Prama RF implants ideal in sinus lift procedures.


Surface

Prima and Prima RF implants have a ZirTi endosseous body (Zirconium Sand-Blasted Acid Etched Titanium), treated with sand-blasting with zirconium oxide and etching with mineral acids. The convergent neck is characterized by a Gold UTM surface (Ultrathin Threaded Microsurface), submitted to a particular controlled passivation process which gives to it a golden yellow colour.

ZirTi Surface - Gold UTM



The UTM neck allows a perfect control of the connection diameter and prevents the plaque accumulation on the connection with the post; moreover, the particular roughness given by the machined neck allows a great adherence of the connection fibres. The colour which characterizes the convergent section allows a natural mimetic of the metal under the soft tissues and under the new materials used in the implantoprosthesis, whose mimetic is due to the translucence and to the transparency.



The ZirTi body is treated with appropriate subtraction techniques, which give the surface the characteristic micromorphology which can remarkably increase the bone-implant contact surface and guarantee a great primary stability.

The validity of the surfaces is documented in numerous **experimental studies**:

Rossi F., Botticelli D., Pantani F., Priscila Pereira F., Salata L.A., Lang N.P.
Bone healing pattern in surgically created circumferential defects around submerged implants: an experimental study in dog
Clin. Oral Impl. Res. 23, 2012, 41-48. doi: 10.1111/j.1600-0501.2011.02170.x

Sivolella S., Bressan E., Salata L.A., Urrutia Z.A., Lang N.P., Botticelli D.
Osteogenesis at implants without primary bone contact – An experimental study in dogs
Clin. Oral Impl. Res. 23, 2012, 542-549 doi: 10.1111/j.1600-0501.2012.02423.x

Rossi F., Lang N.P., De Santis E., Morelli F., Favero G., Botticelli D.
Bone-healing pattern at the surface of titanium implants: an experimental study in the dog
Clin. Oral Impl. Res. 00, 2013, 1-8 doi: 10.1111/clr.12097

Baffone G., Lang N.P., Pantani F., Favero G., Ferri M., Botticelli D.
Hard and soft tissue changes around implants installed in regular-sized and reduced alveolar bony ridges. An experimental study in dogs
Clin. Oral Impl. Res. 00, 2013, 1-6 doi: 10.1111/clr.12306

Key to the implant codes

The implant codes are so-called “mnemonic” codes, i.e. they allow easy identification of the piece. Below is a table showing how the mnemonic codes work using **LA-ZT-425-115** as an example:

type of implant	endosseous morphology	surface	diameter	length
L	A	ZT	425	115
L : Prama implant	A : cylindrical body S : RF tapered body	ZT : ZirTi surface	380 : 3.80 mm 425 : 4.25 mm 500 : 5.00 mm It is the diameter of the implant in its wider point	060 : 6.00 mm 085 : 8.50 mm 100 : 10.00 mm 115 : 11.50 mm 130 : 13.00 mm 150 : 15.00 mm Nominal length which expresses the endosseous length of the implant

All measurements are given in mm, unless indicated otherwise.

Colour codes












A colour code system has been defined in the Prama implant system for identifying the endosseous diameter of the implant (see tables at pages 12 and 13).

The colour code identifies:

- the transfers for the impression taking and the laboratory analogs;
- the final drills;
- the sequence on the surgical tray.










∅ implant	3.80	4.25	5.00
colour code on the pack			

Prama

	3.80	4.25	5.00
6.00	LA-ZT-380-060 	LA-ZT-425-060 	LA-ZT-500-060 
8.50	LA-ZT-380-085 	LA-ZT-425-085 	LA-ZT-500-085 
10.00	LA-ZT-380-100 	LA-ZT-425-100 	LA-ZT-500-100 
11.50	LA-ZT-380-115 	LA-ZT-425-115 	LA-ZT-500-115 
13.00	LA-ZT-380-130 	LA-ZT-425-130 	LA-ZT-500-130 
15.00	LA-ZT-380-150 	LA-ZT-425-150 	LA-ZT-500-150 
vite chirurgica di chiusura	L-VT-340 	L-VT-340 	L-VT-340 

Please note: Nominal length which expresses the endosseous length of the implant. The total length of the implant is 2.80 mm greater than the nominal one, due to the presence of the convergent neck. Each implant is sold with its own surgical cover screw in titanium Gr. 4. The surgical cover screws are also available on sale individually in a sterile pack and must be tightened to 8-10 Ncm.

Prama RF

	3.80	4.25	5.00
6.00	-	LS-ZT-425-060 \varnothing 3.40..... \varnothing 4.25..... 	LS-ZT-500-060 \varnothing 3.40..... \varnothing 5.00..... 
8.50	LS-ZT-380-085 \varnothing 3.40..... \varnothing 3.80..... \varnothing 2.25..... 	LS-ZT-425-085 \varnothing 3.40..... \varnothing 4.25..... \varnothing 2.65..... 	LS-ZT-500-085 \varnothing 3.40..... \varnothing 5.00..... \varnothing 3.40..... 
10.00	LS-ZT-380-100 \varnothing 3.40..... \varnothing 3.80..... \varnothing 2.25..... 	LS-ZT-425-100 \varnothing 3.40..... \varnothing 4.25..... \varnothing 2.65..... 	LS-ZT-500-100 \varnothing 3.40..... \varnothing 5.00..... \varnothing 3.40..... 
11.50	LS-ZT-380-115 \varnothing 3.40..... \varnothing 3.80..... \varnothing 2.25..... 	LS-ZT-425-115 \varnothing 3.40..... \varnothing 4.25..... \varnothing 2.65..... 	LS-ZT-500-115 \varnothing 3.40..... \varnothing 5.00..... \varnothing 3.40..... 
13.00	LS-ZT-380-130 \varnothing 3.40..... \varnothing 3.80..... \varnothing 2.25..... 	LS-ZT-425-130 \varnothing 3.40..... \varnothing 4.25..... \varnothing 2.65..... 	LS-ZT-500-130 \varnothing 3.40..... \varnothing 5.00..... \varnothing 3.40..... 
15.00	LS-ZT-380-150 \varnothing 3.40..... \varnothing 3.80..... \varnothing 2.25..... 	LS-ZT-425-150 \varnothing 3.40..... \varnothing 4.25..... \varnothing 2.65..... 	LS-ZT-500-150 \varnothing 3.40..... \varnothing 5.00..... \varnothing 3.40..... 
vite chirurgica di chiusura	L-VT-340 	L-VT-340 	L-VT-340 

Please note: lNominal length which expresses the endosseous length of the implant. The total length of the implant is 2.80 mm greater than the nominal one, due to the presence of the convergent neck.

Each implant is sold with its own surgical cover screw in titanium Gr. 4. The surgical cover screws are also available on sale individually in a sterile pack and must be tightened to 8-10 Ncm.

See technical characteristics of Gr. 4 titanium on page 13.

Surgical kit

The Prama surgical kit includes all the instruments needed to insert both the standard Prama implants, with cylindrical endosseous morphology, and the Prama RF implants, with tapered body. Prama short implants of cylindrical morphology and height 6.00 mm must be inserted with Drilling Kit Shorty, while Prama short implants of conical morphology and height 6.00 mm must be inserted with Drilling Kit Syra Short (see the following pages). Each type of preparation has the related dedicated drills, whose use sequence is given by coloured marks for the various implant diameters. For the Prama RF there are in the kit also the titanium replies which allow to evaluate the congruity of the receiving site compared to the implant. Together with the kit also templates are supplied, with the graphical representation of the implants, both in real dimension and enlarged of 20% and 30% in order to allow the choice of the implants in their most appropriate dimensions by means of radiographic or tomographic analysis.

The compact dimensions of the kit make it very practical in everyday use and in transport



A practical ratchet is also included that acts as a dynamometric key for the torque of the prosthetic screws and as a surgical key for inserting the implants. The ratchet has a very small head, making it easy to use even in distal sectors

The kit consists of a practical box in Radel with a surgical tray inside that is set-up to hold the instruments according to a guided procedure. The sequences of use of the instruments are indicated by coloured marks

description	code
<p>Complete grommetless surgical kit of the instruments necessary for Prama and Prama RF implants</p>	<p>ZPRAMA-INT</p> 
<p>Radel instrument grommetless tray for Prama and Prama RF instruments</p>	<p>L-TRAY-INT</p> 

Important warning

The surgical kit also contains a test implant (non sterile) which is not to be clinically used, it can be distinguished from the others as it is entirely anodised in blue; it is recommended to use this implant for making trials on the model before starting to use the implants for clinical use, in order to get to know the implant system and its instruments.



1

Precision drill
FS-230



2

Pilot drill
FPT3-200-LXS ◦



3

Intermediate drills
FG-200/280XS ●
FG-330/425XS ●



4

Stop for conical drills
SH-STOP4-FK380 ●
SH-STOP4-FK425 ●
SH-STOP4-FK500 ●

5

Conical drills
SH-FK380-085 ●
SH-FK380-100 ●
SH-FK380-115 ●
SH-FK380-130 ●
SH-FK380-150 ●



27

Stop for pilot drill
STOP4-200-085 ◦
STOP4-200-100 ◦
STOP4-200-115 ◦
STOP4-200-130 ◦
STOP4-200-150 ◦



26

Extension
PROF-CAL3



25

Parallelism pins
PP-2/28



24

Extension
BPM-15



23

Adaptor
AVV-CA-DG-EX



22

Dynamometric ratchet
CRI5-KIT



21

Bone taps
SH-MS-380-CA ●
SH-MS-425-CA ●
SH-MS-500-CA ●




20

Surgical screwdrivers
HSMXS-20-DG HSM-20-DG HSML-20-DG



19

Prosthetic screwdrivers
HSM-20-EX HSML-20-EX HSM-20-CA



6



- Reply**
- SH-380-085-RP ●
 - SH-380-100-RP ●
 - SH-380-115-RP ●
 - SH-380-130-RP ●
 - SH-380-150-RP ●

7



- Conical drills**
- SH-FK425-085 ●
 - SH-FK425-100 ●
 - SH-FK425-115 ●
 - SH-FK425-130 ●
 - SH-FK425-150 ●

8



- Reply**
- SH-425-085-RP ●
 - SH-425-100-RP ●
 - SH-425-115-RP ●
 - SH-425-130-RP ●
 - SH-425-150-RP ●

9



- Conical drills**
- SH-FK500-085 ●
 - SH-FK500-100 ●
 - SH-FK500-115 ●
 - SH-FK500-130 ●
 - SH-FK500-150 ●

10



- Reply**
- SH-500-085-RP ●
 - SH-500-100-RP ●
 - SH-500-115-RP ●
 - SH-500-130-RP ●
 - SH-500-150-RP ●

11



- Final cylindrical drills**
- FFT3-300-LXS ●
 - FFT3-340-LXS ●
 - FFT3-425-LXS ●

12



- Stop for cylindrical drills**
- STOP4-300-085 ●
 - STOP4-300-100 ●
 - STOP4-300-115 ●
 - STOP4-300-130 ●
 - STOP4-300-150 ●

13



- Stop for cylindrical drills**
- STOP4-340-085 ●
 - STOP4-340-100 ●
 - STOP4-340-115 ●
 - STOP4-340-130 ●
 - STOP4-340-150 ●

14



- Stop for cylindrical drills**
- STOP4-425-085 ●
 - STOP4-425-100 ●
 - STOP4-425-115 ●
 - STOP4-425-130 ●
 - STOP4-425-150 ●

15

- Drivers**
- BC-EX230
 - BL-EX230



18

- Prosthetic screwdrivers for Full Head screws**
- L-HSM-EX
 - L-HSML-EX
 - L-HSMXL-EX
 - L-HSM-CA




17

- Easy Insert drivers**
- EASY4-EX230-EX
 - EASYC4-EX230-CA
 - EASYL4-EX230-CA



16

- Bone taps**
- A-MS-380-CA ●
 - A-MS-425-CA ●
 - A-MS-500-CA ●



General indications

The surgical instruments designed for use with the implant systems manufactured by Sweden & Martina are reusable medical devices intended for transient use in the oral cavity for temporary use (not more than 60 minutes), re-usable. The functions of the surgical instruments are to prepare sites for Sweden & Martina implants, to insert the implants in the sites, to tighten and unscrew all the connecting screws (cover screws, healing abutments, screws for posts, abutments, prosthetic screws, transfer screws, etc.).

The surgical instruments manufactured by Sweden & Martina are designed for use with dental implants manufactured by Sweden & Martina. Use of surgical instruments for implant work other than those manufactured by Sweden & Martina limits the responsibility of Sweden & Martina and renders the product warranty void. Sweden & Martina declines all responsibility for use of any non-original instruments. Sweden & Martina surgical instruments are sold in NON-STERILE packs. Before use, they must be cleaned, disinfected and sterilised according to the instructions reported below. Failure to follow these warnings may expose the patient to infection. The materials used for manufacturing the surgical instruments manufactured by Sweden & Martina were selected based on the properties indicated for their intended use according to directive 93/42, implemented in Italy with Law 46/97, Annex I – Essential Requirements, point 7.1.

Each packaging indicates the code, description of the contents and batch number. These same details, which are also indicated on the labels inside the packs, must always be provided by the practitioner in any relevant correspondence.

All the devices are identified by an instrument code, which is laser marked onto the body of each instrument. If there is not enough space to include the full code, the elements for unequivocally identifying the device (e.g. diameter or length) are provided. When handling the devices, both during use and during cleaning and sterilisation, it is recommended to use surgical gloves for personal protection from bacterial contaminations. Failure to follow these instructions may cause cross-infection.

Key to the implant codes: surgical instruments

The implant codes are so-called “mnemonic” codes, i.e. they allow easy identification of the piece. Below is a table showing how the mnemonic codes work using different types of instruments as an example.

examples	type of component	diameter	length
The range of instruments is vast, we indicate some examples of the main families of instruments.	-	Normally it is the \varnothing of the implant for the insertion or of the preparation for which the instrument should be used.	This measurement is normally linked to the height of the component, or to other important measurements that characterise it, or it is a letter which defines whether a post is repositionable or not.
FFT3-300-LXS	FFT3 : final cylindrical drill	300: 3.00 mm, for the preparation of the \varnothing 3.80 mm implant	115: 11.50 mm
STOP4-200-085	STOP4 : stop for cylindrical drills	200: 2.00 mm	085: 0.85 mm
PP-2/28	PP : parallelism pin	2/28: from 2.00 mm to 2.80 mm	-

Drills

All Sweden & Martina drills are made of **stainless steel** with **high resistance to corrosion and wear**. They are intended for mechanical use, i.e. they have a shank with a right angle attachment and must be used with a suitable micromotor. The extreme accuracy of design and production allows use completely free from vibrations and oscillations. However, incorrect insertion of the instruments in the handpiece will cause instrument vibration, eccentric rotation, early wear and shaft buckling. Suitable surgical micromotors only should be used. Micromotors should be checked regularly by their manufacturers, according to the indications given by the manufacturers, to prevent potential malfunctions (e.g. axle shifts for transmission shafts, worn or faulty forceps, etc.).

Failure to follow the instructions provided may cause surgical complications and consequent damage to the patient's health. It is recommended to use the rotation speeds indicated in the procedures on page 50 to prevent the development of bone necrosis. Lever movements increase the risk of instrument breakage and should therefore be avoided. Changes in speed should be avoided in general. Never apply pressure such as to force the instrument to stop rotating. This could lead to an excessive increase in heat in the tissues being drilled, with consequent bone necrosis, and damage both the instrument and the appliance (micromotor) used. This could also lead to breakage of the instrument. Using an intermittent approach, with a back and forth movement in a vertical direction, prevents overheating and wear of the working part and an undesirable increase in the temperature in the tissues being cut.

Suitable coolant should be used. Inadequate irrigation can lead to bone necrosis.

Drill wear depends to a large extent on the type and density of the drilled bone: harder bone leads to greater instrument wear. For greater safety and caution, given the device's capacity for resistance to wear, drills should not be used for more than **20 work cycles** and should be replaced earlier if the instruments lose their cutting ability. These recommended 20 cycles should be considered a rough guide. Always check the instrument's residual cutting capacity after each procedure. Sweden & Martina decline responsibility for the use of blunted instruments.

Never sharpen drills before use. Never use damaged, buckled or worn instruments.



Precision drill FS-230

The Precision drill is made of surgical stainless steel. It is used to cut the cortical bone, so it is very sharp and pointed. The design of the blades ensures efficient cutting with both the tip and the edge. It has a maximum diameter of 2.30 mm. The laser marking at 4.80 mm indicates the depth to which the drill should always be inserted to obtain a suitable guiding hole for the next drills.



Important warning

The Precision drill comes with a protective silicone sheath. The sole purpose of this protective sheath is to protect the instrument during transportation and it must be removed before first use. Since this drill is extremely sharp, special caution is required during handling.


Pilot drill FPT3-200-LXS

The pilot drill \varnothing 2.00 is used to prepare the initial hole for preparing the site. The drill is easy to identify, thanks to the presence of a white ring and to the code laser-etched on the drill shank. It has laser-etched depth marks, a cylindrical shape and a spiral with two cutting edges. It must be used with abundant external irrigation.



LL: Total length of the working part, including the tip.

LS: Length of the tip. This measurement must be calculated in addition to the length of the preparation hole.



Important warning

The drills always make a hole that is longer than the implant to be inserted. The oversizing (LS) is equal to the height of the tip of the drill that is being used.

code	∅	LS	LL
FPT3-200-LXS	2.00	0.58	19.30

Pilot drill stops

Stops are devices to be fitted in tip → shank direction on drills suited to receive them. They make it possible to restrict the working length of a drill to a pre-set height.


height	8.50 mm	10.00 mm	11.50 mm	13.00 mm	15.00 mm
stop	STOP4-200-085	STOP4-200-100	STOP4-200-115	STOP4-200-130	STOP4-200-150
					

Always check that the stop is inserted at the desired height. Incomplete insertion may reduce the preparation height. Any insertion difficulties can be resolved by loosening the stop tabs slightly, using forceps. It is also recommended to check the retention exerted by the stop, as if retention is too weak the instrument will fall off the drill during operation. In the event of reduced retention capacity, simply tighten the tabs by hand or using forceps.



Intermediate drills

Intermediate drills are drills with two cutting edges suitable for progressively widening the preparations in relation to the diameter of the drills to be used in succession. They have two small steps with an initial guide with a progressive diameter and final diameter, respectively equal to 2.00/2.80 and 3.30/4.25 mm. They have reference laser markings that range from a height of 8.50 to 10.00 mm. For shorter preparations, they must be used until the end stop (the guide is not a cutting edge).




Other intermediate drills

Optionally two \varnothing 2.50 mm stainless steel intermediate drills are available. These drills, in cylindrical version for Prama implants and conical version for Prama RF implants, are useful for underpreparation protocols. The related stops are available, meant to guarantee a safe preparation.




Both the drills are not contained in the Prama surgical kit. The \varnothing 2.50 mm cylindrical drill and the related stops are contained in the kit on sale with the code **KIT-INTEGRA-F250**, while the \varnothing 2.50 mm conical drill and the related stops are contained in the kit on sale with the code **KIT-INTEGRA-SE-FK250**. Drills and stops are also available separately.

Please note: the drills always make a hole that is longer than the implant to be inserted. The oversizing (LS) is equal to the height of the tip of the drill that is being used. See drawing on the side.



LL: Total length of the working part, including the tip.

LS: Length of the tip. This measurement must be calculated in addition to the length of the preparation hole.



cylindrical drill ϕ 2.50 mm and related stops

<p>FFT3-250-LXS cylindrical drill</p> <p>ϕ 2.50 19.50 0.72</p>	<p>STOP4-250-070 stop 7.00 mm for cylindrical drill</p>	<p>STOP4-250-085 stop 8.50 mm for cylindrical drill</p>	<p>STOP4-250-100 stop 10.00 mm for cylindrical drill</p>	<p>STOP4-250-115 stop 11.50 mm for cylindrical drill</p>	<p>STOP4-250-130 stop 13.00 mm for cylindrical drill</p>	<p>STOP4-250-150 stop 15.00 mm for cylindrical drill</p>
---	--	--	---	---	---	---

conical drill ϕ 2.50 mm and related stops*


<p>SE-FK250 conical drill</p> <p>ϕ 2.50 18.00 ϕ 1.80 0.50</p>	<p>CSR-STOP-2028-085 stop 8.50 mm for conical drill</p>	<p>CSR-STOP-2028-100 stop 10.00 mm for conical drill</p>	<p>CSR-STOP-2028-115 stop 11.50 mm for conical drill</p>	<p>CSR-STOP-2028-130 stop 13.00 mm for conical drill</p>	<p>CSR-STOP-2028-150 stop 15.00 mm for conical drill</p>
--	--	---	---	---	---

Final cylindrical drills and related stops

Made of stainless steel with high resistance to corrosion and wear, the Prama final drills have a number of cutting edges proportioned to the diameter of the hole, in order to allow a continue and homogeneous cutting movement and a better stability of the instrument during the surgical phases. All this allows to obtain high precision implant preparations, with following easiness in the insertion of the implant.

The use of these drills is recommended with the related depth stops, also included in the kit.

LL: Total length of the working part, including the tip.
LS: Length of the tip. This measurement must be calculated in addition to the length of the preparation hole.



∅ implant	3.80 mm	4.25 mm	5.00 mm
cylindrical drills	FFT3-300-LXS	FFT3-340-LXS	FFT3-425-LXS
	 19.60 0.87	 19.70 0.98	 20.00 1.23

Important warning

The drills always make a hole that is longer than the implant to be inserted. The oversizing (LS) is equal to the height of the tip of the drill that is being used.


Stop per frese cilindriche

\varnothing implant	3.80 mm	4.25 mm	5.00 mm
stop for preparations H 8.50 mm	STOP4-300-085 	STOP4-340-085 	STOP4-425-085 
stop for preparations H 10.00 mm	STOP4-300-100 	STOP4-340-100 	STOP4-425-100 
stop for preparations H 11.50 mm	STOP4-300-115 	STOP4-340-115 	STOP4-425-115 
stop for preparations H 13.00 mm	STOP4-300-130 	STOP4-340-130 	STOP4-425-130 
stop for preparations H 15.00 mm	STOP4-300-150 	STOP4-340-150 	STOP4-425-150 

It is recommended always to check that the stop is inserted at the desired height. Incomplete insertion may reduce the preparation height. Any insertion difficulties can be resolved by loosening the stop tabs slightly, using forceps. It is also recommended to check the retention exerted by the stop, as if retention is too weak the instrument will fall off the drill during operation. In the event of reduced retention capacity, simply tighten the tabs by hand or using forceps.

Final conical drills and related stops

The conical drills for Prama RF implants are also made of stainless steel with high resistance to corrosion and wear. They present a number of cutting edges proportional to the hole diameter, so as to allow a continuous and homogeneous cutting movement and greater instrument stability during operation. They are distinguished by a coloured ring that makes it easy to recognize the instruments intended for each diameter. All this results in very precise implant preparations, which are the key to success of conical implants. The kit contains 15 conical drills, each one of which forms the final hole for the implant with diameter and height referred to by the instrument code.




LT: Total length of the working part, including the shank

LL: Total length of the working part, including the tip

LS: Length of the overpreparation

Important warning

The drills always make a hole that is longer than the implant to be inserted. The oversizing (LS) is equal to the difference between the length of the working part of the drill and the nominal height of the implant. For details of the sizes of the different drills, refer to the table on the next page.

\varnothing implant	3.80 mm	4.25 mm	5.00 mm
8.50	<p>SH-FK380-085</p>  <p>\varnothing 3.60 \varnothing 2.20 8.50 0.42</p>	<p>SH-FK425-085</p>  <p>\varnothing 4.00 \varnothing 2.60 8.50 0.54</p>	<p>SH-FK500-085</p>  <p>\varnothing 4.75 \varnothing 3.35 8.50 0.75</p>
10.00	<p>SH-FK380-100</p>  <p>\varnothing 3.60 \varnothing 2.20 10.00 0.44</p>	<p>SH-FK425-100</p>  <p>\varnothing 4.00 \varnothing 2.60 10.00 0.56</p>	<p>SH-FK500-100</p>  <p>\varnothing 4.75 \varnothing 3.35 10.00 0.77</p>
11.50	<p>SH-FK380-115</p>  <p>\varnothing 3.60 \varnothing 2.20 11.50 0.46</p>	<p>SH-FK425-115</p>  <p>\varnothing 4.00 \varnothing 2.60 11.50 0.57</p>	<p>SH-FK500-115</p>  <p>\varnothing 4.00 \varnothing 2.60 11.50 0.79</p>
13.00	<p>SH-FK380-130</p>  <p>\varnothing 3.60 \varnothing 2.20 13.00 0.47</p>	<p>SH-FK425-130</p>  <p>\varnothing 4.00 \varnothing 2.60 13.00 0.59</p>	<p>SH-FK500-130</p>  <p>\varnothing 4.75 \varnothing 3.35 13.00 0.80</p>
15.00	<p>SH-FK380-150</p>  <p>\varnothing 3.60 \varnothing 2.20 15.00 0.52</p>	<p>SH-FK425-150</p>  <p>\varnothing 4.00 \varnothing 2.60 15.00 0.64</p>	<p>SH-FK500-150</p>  <p>\varnothing 4.75 \varnothing 3.35 15.00 0.85</p>

SURGICAL INSTRUMENTS

drill code	corresponding implant	nominal \emptyset	minimum \emptyset	maximum \emptyset	LT	LL	LS	colour code
SH-FK380-085	LS-ZT-380-085	3.80	2.20	3.60	30.92	8.92	0.42	verde
SH-FK380-100	LS-ZT-380-100	3.80	2.20	3.60	32.44	10.44	0.44	verde
SH-FK380-115	LS-ZT-380-115	3.80	2.20	3.60	33.96	11.96	0.46	verde
SH-FK380-130	LS-ZT-380-130	3.80	2.20	3.60	35.47	13.47	0.47	verde
SH-FK380-150	LS-ZT-380-150	3.80	2.20	3.60	37.52	15.52	0.52	verde
SH-FK425-085	LS-ZT-425-085	4.25	2.60	4.00	31.04	9.04	0.44	blu
SH-FK425-100	LS-ZT-425-100	4.25	2.60	4.00	32.56	10.56	0.56	blu
SH-FK425-115	LS-ZT-425-115	4.25	2.60	4.00	34.07	12.07	0.57	blu
SH-FK425-130	LS-ZT-425-130	4.25	2.60	4.00	35.59	13.59	0.59	blu
SH-FK425-150	LS-ZT-425-150	4.25	2.60	4.00	37.64	15.64	0.64	blu
SH-FK500-085	LS-ZT-500-085	5.00	3.35	4.75	31.26	9.25	0.75	magenta
SH-FK500-100	LS-ZT-500-100	5.00	3.35	4.75	32.77	10.77	0.77	magenta
SH-FK500-115	LS-ZT-500-115	5.00	3.35	4.75	34.29	12.29	0.79	magenta
SH-FK500-130	LS-ZT-500-130	5.00	3.35	4.75	35.80	13.80	0.80	magenta
SH-FK500-150	LS-ZT-500-150	5.00	3.35	4.75	37.85	15.85	0.85	magenta

Stops for conical drills

The kit contains a stop for each diameter of the final conical drills, for **inserting the drill from the tip**. They are suitable for limiting the working length to predetermined heights. With the same working diameter, the same stop is compatible with all the drill lengths, as explained in the following table:

stop for conical drills	SH-STOP4-FK380	SH-STOP4-FK425	SH-STOP4-FK500
			
colour code	green	blue	magenta
nominal \varnothing corresponds to the implant diameter	3.80	4.25	5.00
drill for implant H 8.50 mm	SH-FK380-085	SH-FK425-085	SH-FK500-085
drill for implant H 10.00 mm	SH-FK380-100	SH-FK425-100	SH-FK500-100
drill for implant H 11.50 mm	SH-FK380-115	SH-FK425-115	SH-FK500-115
drill for implant H 13.00 mm	SH-FK380-130	SH-FK425-130	SH-FK500-130
drill for implant H 15.00 mm	SH-FK380-150	SH-FK425-150	SH-FK500-150

As already indicated with regard to the pilot drill stops, in this case too it is recommended always to check that the stop is inserted at the desired height. Incomplete insertion may reduce the preparation height. Any insertion difficulties can be resolved by loosening the stop tabs slightly, using forceps. It is also recommended to check the retention exerted by the stop, as if retention is too weak the instrument will fall off the drill during operation. In the event of reduced retention capacity, simply tighten the tabs by hand or using forceps.

Reply: replies for Prama RF implants

The Reply replies are made of Gr. 5 titanium and reply the morphology of the final drills of the related Prama RF conical implants. They are useful to verify the depth of the preparation hole made with the final drills, and to verify the axis of the preparation made with the drill.





\varnothing implant	3.80 mm	4.25 mm	5.00 mm
H 8.50 mm	<p>SH-380-085-RP</p>  <p>\varnothing 3.60 \varnothing 2.30 8.50</p>	<p>SH-425-085-RP</p>  <p>\varnothing 4.00 \varnothing 2.70 8.50</p>	<p>SH-500-085-RP</p>  <p>\varnothing 4.75 \varnothing 3.45 8.50</p>
H 10.00 mm	<p>SH-380-100-RP</p>  <p>\varnothing 3.60 \varnothing 2.30 10.00</p>	<p>SH-425-100-RP</p>  <p>\varnothing 4.00 \varnothing 2.70 10.00</p>	<p>SH-500-100-RP</p>  <p>\varnothing 4.75 \varnothing 3.45 10.00</p>
H 11.50 mm	<p>SH-380-115-RP</p>  <p>\varnothing 3.60 \varnothing 2.30 11.50</p>	<p>SH-425-115-RP</p>  <p>\varnothing 4.00 \varnothing 2.70 11.50</p>	<p>SH-500-115-RP</p>  <p>\varnothing 4.75 \varnothing 3.45 11.50</p>
H 13.00 mm	<p>SH-380-130-RP</p>  <p>\varnothing 3.60 \varnothing 2.30 13.00</p>	<p>SH-425-130-RP</p>  <p>\varnothing 4.00 \varnothing 2.70 13.00</p>	<p>SH-500-130-RP</p>  <p>\varnothing 4.75 \varnothing 3.45 13.00</p>
H 15.00 mm	<p>SH-380-150-RP</p>  <p>\varnothing 3.60 \varnothing 2.30 15.00</p>	<p>SH-425-150-RP</p>  <p>\varnothing 4.00 \varnothing 2.70 15.00</p>	<p>SH-500-150-RP</p>  <p>\varnothing 4.75 \varnothing 3.45 15.00</p>

Countersink drills


Countersink drills are optionally available. They allow to prepare the seat for the most coronal part of Prama and Prama RF implants in presence of a very thick cortical bone. The drills are characterized by an apical portion that is guided into the hole created by the final cylindrical drill and by a standard height of the working part of 6.60 mm.



Drills are available with 3.80, 4.25 and 5.00 mm diameters

∅ implant	3.80 mm	4.25 mm	5.00 mm
countersink drills	A-FCC-380 	A-FCC-425 	A-FCC-500 


The particular morphology of the countersink drills allows to widen the initial part of the hole bored by the final cylindrical drills to prepare the cortical bone portion, that will be in contact with the neck of the implant. The maximum recommended speed is 1,000 rpm. Each drill has to be used only for the implant of equal diameter.



Drills for distal sectors

As an option, shorter drills are available that are very practical in distal sectors with limited oral opening. They come in a wide range of diameters and are also useful for preparations in extremely compact bone where, in the most coronal portion, it is desired to widen the preparation diameter by 0.10 mm with respect to the size of the standard drills to facilitate the insertion of the implants. On the other hand, in low-density bone they can be used to under-prepare the implant site so as to obtain optimum primary stability.

- LL:** Total length of the working part, including the tip.
- LS:** Length of the tip. This measurement must be calculated in addition to the length of the preparation hole.




Important note: the drills always make a hole that is longer than the implant to be inserted. The oversizing (Ls) is equal to the height of the tip of the drill that is being used. See drawing on the side.

∅ drill	2.00	2.80	2.90	3.00	3.20
drills for distal sector	FPT5-200-LXS	FFT5-280-LXS	FFT5-290-LXS	FFT5-300-LXS	FFT5-320-LXS
∅ frese	3.30	3.40	3.60	4.25	4.45
drills for distal sector	FFT5-330-LXS	FFT5-340-LXS	FFT5-360-LXS	FFT5-425-LXS	FFT5-445-LXS

The drills for distal sectors are not included in any surgical kit and must be purchased separately. They cannot be used with depth stops.

Bone taps

The Prama and Prama RF implants are self-tapping implants with excellent cutting and insertion capabilities; however, the use of a bone tap is recommended in all cases where the type of bone requires it in order to facilitate the insertion of the fixture. They are available both with right angle shank and with a connector for dynamometric ratchet and have a specific design with two different endosseous morphology.



∅ implant	3.80 mm	4.25 mm	5.00 mm
Bone taps with right angle attachment for Prama implant	A-MS-380-CA  ∅ 2.30	A-MS-425-CA  ∅ 2.50	A-MS-500-CA  ∅ 3.55
Bone taps with connection for ratchet for Prama implants*	A-MS-380  ∅ 2.30	A-MS-425  ∅ 2.50	A-MS-500  ∅ 3.55
Short bone taps with connection for ratchet for Prama implants*	A-MSC-380  ∅ 2.25	A-MSC-425  ∅ 2.45	A-MSC-500  ∅ 3.50
Bone taps with right angle attachment for Prama RF implants	SH-MS-380-CA  ∅ 2.55	SH-MS-425-CA  ∅ 2.65	SH-MS-500-CA  ∅ 3.40
Bone taps with connection for ratchet for Prama RF implants*	SH-MS-380  ∅ 2.55	SH-MS-425  ∅ 2.65	SH-MS-500  ∅ 3.40

*Optional instruments not included in the surgical kit, to be ordered separately.

Osteotomes


Osteotomes for expansion protocols are available for each Prama and Prama RF endosseous morphology, not included in the surgical kit. The laser-etched codes on the handles indicate the osteotome diameter, to make it easier to recognize the correct surgical sequence. A practical universal instrument case for storing and organising them is available and can be ordered separately using the code **OS-TRAY**.

Prama osteotomes

Prama osteotomes present laser markings of all the available heights.

code	E-OS-020-PP	E-OS-090-PP	E-OS-160-PC	E-OS-200-PC	E-OS-240-PC
					
 ϕ 0.20 ϕ 0.90 ϕ 1.60 ϕ 2.00 ϕ 2.40
description	Osteotome ϕ 0.20 flat tip	Osteotome ϕ 0.90 flat tip	Osteotome ϕ 1.60 concave tip	Osteotome ϕ 2.00 concave tip	Osteotome ϕ 2.40 concave tip

Osteotomes are optional instruments that are not included in the surgical kit. They can be purchased separately and singularly.




















description	code
Universal case in Radel for osteotomes for Prama and Prama RF Can hold up to 12 instruments	OS-TRAY



Prama RF osteotomes

Prama RF osteotomes are designed considering the height and diameter of the implant to be inserted. They present two laser-etched markings: one corresponds to the height of the fixture and the other one, at a lower height, is useful for an intermediate check during the preparation phase. In the osteotome for 10.00 mm implants the laser markings correspond to the height 8.50 and 10.00 mm, therefore the instrument can also be used for inserting implants of height 8.50 mm.

ϕ implant	3.80		4.25		5.00	
osteotome for implants H 8.50 and 10.00 mm	SH-OS-380-100-PP 	SH-OS-380-100-PR 	SH-OS-425-100-PP 	SH-OS-425-100-PR 	SH-OS-500-100-PP 	SH-OS-500-100-PR 
osteotome for implants H 11.50 mm	SH-OS-380-115-PP 	SH-OS-380-115-PR 	SH-OS-425-115-PP 	SH-OS-425-115-PR 	SH-OS-500-115-PP 	SH-OS-500-115-PR 
osteotome for implants H 13.00 mm	SH-OS-380-130-PP 	SH-OS-380-130-PR 	SH-OS-425-130-PP 	SH-OS-425-130-PR 	SH-OS-500-130-PP 	SH-OS-500-130-PR 
osteotome for implants H 15.00 mm	SH-OS-380-150-PP 	SH-OS-380-150-PR 	SH-OS-425-150-PP 	SH-OS-425-150-PR 	SH-OS-500-150-PP 	SH-OS-500-150-PR 
tip	flat	round	flat	round	flat	round

Osteotomes are optional instruments that are not included in the surgical kit. They can be purchased separately and singularly.

Parallelism pins


Parallelism pins can be used to check the insertion axis of the implants and the parallelism between more fixtures. They have one side with diameter \varnothing 2.00 mm and the other \varnothing 2.80 mm, so as to be used after the drills with the same diameter.



description	code
parallelism pin with a \varnothing 2.00 mm side and a \varnothing 2.80 mm	PP-2/28 

Parallelism pins with depth lines are available optionally, they allow the control of the preparation depth during the first surgical step, thanks to the presence of dedicated lines in the side with \varnothing 2.00 mm. As the lines have a reduced diameter in comparison with the pin body, it is possible to distinguish them also on the x-ray images. The other side of the instruments has a diameter of \varnothing 2.80 mm and presents a hole for safety thread.

The small version of the pin, the one with a shorter \varnothing 2.80 side, is useful for patients with limited oral opening or the use in distal sectors.



description	code
parallelism pin with depth lines, large	PPTL-2-28 
parallelism pin with depth lines, small	PPTS-2-28 


Important warning

It is recommended to pass a thread through the hole in the centre of the pin to prevent it falling.

Parallelism pins with depth lines are not contained in the surgical kit, they can be ordered singularly for separate.




Depth gauge PROF3

It is a practical instrument which allows to verify the depth of the preparation holes and the distance between the implants. It is not included in any surgical kit, it can be ordered separately.



X-ray templates

The surgical kits also contain templates for the graphic representation of the implant measurements to allow choosing the most suitable implant diameters and lengths by means of radiographic or tomographic methods. The templates are available in three versions: in real dimension and enlarged of 20% and 30% .

description	code
X-ray template for Prama and Prama RF implants, real dimensions	<p>L-L100</p> 
X-ray template for Prama and Prama RF implants, dimensions increased by 20%	<p>L-L120</p> 
X-ray template for Prama and Prama RF implants, dimensions increased by 30%	<p>L-L130</p> 

Easy Insert driver

The Prama implant does not require a mounter for inserting into the implant site because it is engaged directly inside the connection by practical Easy Insert drivers designed to guarantee a safe grip, to prevent deformations to the corners of the connections and at the same time to allow easy removal from the implant wells. The use of these drivers makes the surgical procedure of insertion extremely easy.


The Easy Insert drivers aid **visibility of the operating field**, do not occupy much space, and allow adjusting the connection hexagon properly because they present a hexagonal visual index corresponding to the prosthetic index and three laser marks on three of the six sides

The presence of a wide hexagon allows **engaging the patented Easy Insert drivers easily and safely** for insertion of the implants into the relative sites




A single instrument that allows the insertion of all Prama implant diameters

The whole is extremely safe and reliable with the use of a special **titanium O-ring** that engages inside the connection



The photo on the left shows how a traditional instrument (in green) edges itself into the connection (in grey). This geometry inevitably determines the grip and deformation of the actual session



The **special design** of the Easy Insert drivers prevent any deformations of the implant connection, since the driver's faces (not the corners) are the ones in contact with the the connection walls

The **dodecagonal design** of the drivers prevents deformations to the implant connection, thus guaranteeing extremely high prosthetic stability and precision.

When using the Easy Insert with ratchet, as when using any other instrument for inserting implants with a dynamometric key, it is likewise advisable to take care to keep the working axis as perpendicular as possible.

It is also fundamental for the movement performed with the ratchet during tightening to be slow and uniform, avoiding brusque movements as much as possible. If these recommendations are not respected and the insertion torque is exceeded, the instrument should get broken: for this reason a preferential failure point is present above the black laser marks, to help the operator to easily remove the driver.

It is recommended to grip the ratchet in the part closest to the connection and to maintain a light and constant pressure with one finger, to allow greater stability during tightening.

description	code	included in the kit
short driver with right angle shank	EASYC4-EX230-CA 	Prama
long driver with right angle shank	EASYL4-EX230-CA 	Prama
driver with connector for dynamometric key	EASY4-EX230-EX 	Prama

Maintenance and care of the Easy Insert drivers

The Easy Insert drivers are supplied pre-mounted with the special titanium O-rings. Since they are mechanical components, the retainer rings are subject to wear over time and can lose their elasticity and functionality.

The O-rings cannot be replaced, but it is necessary to replace the instrument. The Easy Inserts were tested to be good for 40 uses in the worst conditions of use. These limits can therefore change depending on the conditions of use.

However, it is always a good idea to check its good functionality even during the cleaning and sterilisation operations. For this reason and to allow the doctor to familiarise himself with the Easy Inserts, the surgical kit contains a “test implant” that has not been treated or sterilised; it can be distinguished from the others as it is in blue.

Important warning

It is recommended to use the Easy Insert with a torque between 50 Ncm and 70 Ncm. Thanks to tests performed on models, it has been observed that from 70 Ncm to 100 Ncm slight frictions between the instrument and the implant connection are possible, they can be avoided with a slight shaking movement of the Easy Insert in the connection. From 100 Ncm to 200 Ncm higher frictions are possible, they can be solved with a simple counter-rotation movement (at 40 Ncm) in order to remove the instrument from the connection. It is moreover recommended to end the bone tapping phase using a dynamometric ratchet.

Drivers for fixation screws

All the drivers are made of stainless steel for surgical use. There are two types of drivers for the Prama implants: the traditional ones (on the left in the picture) and the ones for the screws with Full Head technology (on the right). They differ in the design of the tip, studied in the first case to join a screw with internal hexagonal connection and in the other with external hexagonal connection, therefore they are not interchangeable. In both cases the slightly conical coupling between the driver and the screw allows an appropriate retention when carrying the screw in the oral cavity. Both drivers families are available in different shank lengths, in order to facilitate the ergonomics depending on the patient anatomy. The standard drivers are available also in the hand one-piece version, this means they are integral with the hand knob which allows the grip.

Regularly verify that this functionality have not been lost due to wear.



Important warning

Excessive torques can damage the thread of the well or of the sharp edges of the connecting screws and damage the thread of the drivers, causing also severe intra-surgical and prosthetic complications. The recommended torque for the tightening of the different components are summarized in the following chart:

description	recommended torque
surgical cover screws, healing abutments	(manually) 8-10 Ncm
all the prosthetic screws	20-25 Ncm
all the prosthetic components with direct screwing on the implant	25-30 Ncm
transfer fixation screws	(manually) 8-10 Ncm

Given the importance of the tightening torque, it is recommended to use always the drivers with hexagonal connection, keeping always the exerted torque under control with the ratchet. To facilitate the joint of the screws or of the threaded sections of the prosthetic components, the screwing can be started with the hand drivers.



Important warning

It is recommended to pass a thread through the hole on the top of the knob to prevent it falling.




Surgical screwdrivers

Their design makes them very practical in the surgical phases for the screwing of the surgical connecting screws and for the phases of uncovering and management of the healing abutments. They must not be used in the final prosthetic screws since they do not allow the control tightening torque.

description	code	included in the kit
screwdriver for tap screws and fixation screws, digital, extra-short	HSMXS-20-DG 	Prama
screwdriver for tap screws and fixation screws, digital, short	HSM-20-DG 	Prama
screwdriver for tap screws and fixation screws, digital, long	HSML-20-DG 	Prama





Prosthetic screwdrivers (for screws with Full Head technology)

The specific screwdrivers for the Full Head technique screws are available in the version with hexagonal connection for ratchet, with different shank lengths. A screwdriver with right angle connector is also available. The instrument's tip has a hexagonal notch, which connects the full hexagon of the Full Head screws, giving the retention needed for the carriage of the screw. These screwdrivers cannot be used for the connecting screws of the temporary posts or for other types of screws in the catalogue.

description	code	included in the kit
screwdriver for screws with Full Head technology, digital, extra short	L-HSM-EX 	Prama
screwdriver for screws with Full Head technology, digital, short	L-HSML-EX 	Prama
screwdriver for screws with Full Head technology, digital, long	L-HSMXL-EX 	Prama
screwdriver for screws with Full Head technology, for right angle	L-HSM-CA 	Prama

Prosthetic screwdrivers for standard screws (surgical cover screws, healing abutments, standard prosthetic screws)

The screwdrivers with upper hexagonal connection have been designed to be used with the dynamometric ratchet with function of control of the torque. In the kit there are two versions available the long one and the short one, an optional extra-long version is also available, necessary when the length of the hole for the screw to pass inside the posts is greater than 13.00 mm.

description	code	included in the kit
screwdriver for fixation screws, with connector for dynamometric ratchet or digital connector, short	HSM-20-EX 	Prama
screwdriver for fixation screws, with connector for dynamometric ratchet or digital connector, long	HSML-20-EX 	Prama
screwdriver for fixation screws, with connector for dynamometric ratchet or digital connector, extra-long	HSMXL-20-EX 	Not included in the surgical kit, available separately
screwdriver for fixation screws, with right angle shank	HSM-20-CA 	Prama
screwdriver for fixation screws, with right angle shank, extra-short	HSMXS-20-CA 	Not included in the surgical kit, available separately

Important warning

All the ratchet drivers have a red polymer O-ring in the connecting hexagon that guarantees friction between the instruments and therefore a correct grip of the components. This O-ring must be checked periodically and replaced when worn or when no longer able to exert the correct friction. A kit of 5 spare O-rings is available, which can be ordered with code ORING180-088.

Drivers for intraoperative removal of implants

description	code	included in the kit
short driver	BC-EX230 	Prama
long driver	BL-EX230 	Prama

Other instruments

description	code	included in the kit
Extension for bone taps, mounters, drivers and manual drivers, with hexagonal connector for dynamometric key	BPM-15 	Prama
Extension for surgical drills	PROF-CAL3 	Prama
Mechanical adapter with right angle shank for instruments with hexagonal connector	B-AVV-CA3 	Prama
Driver for right angle and manual instruments and instruments with hexagonal connection for ratchet	AVV-CA-DG-EX 	Not included in the surgical kit, available separately
Digital knob for bone taps, mounters, screwdrivers and drivers	AVV3-MAN-DG 	Prama

Spare O-ring


description	code	included in the kit
Box with 5 spare o-rings for all accessories with hexagonal connector for dynamometric key	ORING180-088 	Not included in the surgical kit, available separately

Mounter

description	code	included in the kit
Mounter	MOU-EX230 	Not included in the surgical kit, available separately
Mounter stop key	CM2 	Not included in the surgical kit, available separately

Dynamometric ratchet CRI5-KIT

The surgical kit of the implant system contains a special ratchet (CRI5-KIT), with its own adjustment key, for quickly screwing the torque adjustment ring nut, and with gel lubricant for maintenance. The ratchet may be used with torque adjustment from 10 to 70 Ncm or in a blocked position without torque control. When using as a prosthetic ratchet for fastening the screws, refer to the torque values given in the table on the previous page. The ratchet key CRI5 is a multi-purpose instrument that can be disassembled, and is sold unsterile.



Before each use, this instrument must be cleaned and sterilised according to the instructions on page 50. Adequate maintenance, performed following in detail all the step by step instructions for the disassembly and correct reassembly of the device during cleaning operations, is essential for the correct functioning of the device and for its durability. Personnel who use this tool must be suitably trained, and they must have read the instructions in this manual prior to handling the device. After sterilisation, the key is ready for use. A test to verify the correct assembly and functioning of the key is necessary before any surgical or prosthetic interventions.

The torque is adjusted by aligning the marking of the desired torque in the circular opening of the handle. The "IN" arrow legible on the top of the head indicates the screwing position of the key. The "OUT" arrow legible on the top of the head indicates the loosening or unscrewing position. An unlimited torque position is obtained by positioning the torque adjustment device up to the line marked "R" on the handle of the ratchet body.



The ring nut may be screwed and unscrewed by hand, but to speed up these operations the kit also contains a driver that allows it to be turned quickly. Any deterioration of the screwing, insertion and torque mechanisms must be checked by personnel responsible for the use and maintenance of this dental instrument. The pieces of this mechanism are not interchangeable; one piece from one key cannot be replaced by a piece from another key as each ratchet is calibrated INDIVIDUALLY. If a piece is lost, please return the instrument to Sweden & Martina for repair. No components for assembling the ratchet can be sold individually. Failure to follow the instructions provided may cause problems of maintenance and stability of the prosthesis.



Important warning

The torque is adjusted by screwing/unscrewing the ring nut located at the bottom of the instrument's handle. The torque must always be adjusted on the rise, starting screwing from a lower value until the desired torque is reached, or unscrewing the ring nut in a clockwise direction. To do this, if it is necessary to set a torque lower than the last one used, you must unscrew the ring nut by two turns below the value of the desired new torque, and work up to that value by rescrewing the ring nut in a clockwise direction.




In order to set a torque value, turn the ring nut in a clockwise direction until the wanted value.

To turn down a torque value of work of the ratchet, first it is necessary to unscrew the ring nut in a anticlockwise direction until reaching a value inferior of the wanted one, then proceed with the clockwise direction screwing until the chosen torque.

Torque wrench with control lever TWL

The surgical kit of the CSR-DAT implant system includes a special dynamometric key with a control lever (TWL). The torque wrench can be used to indicate the value of the torque applied during the surgical phases of screwing and unscrewing, showing values from 10 to 90 Ncm. It is supplied complete with a specific adaptor that allows it to be used with surgical instruments with a hexagonal connection.

The torque wrench with control lever TWL is a multipurpose instrument that can be dismantled, and it is sold as non-sterile.



Every time this instrument is used, it must first be cleaned and sterilized following the instructions on page 52.

Adequate maintenance, carried out by scrupulously following all the steps indicated for dismantling and reassembling of the torque wrench during cleaning operations is essential for its correct use and to prolong its shelf life. The personnel using this instrument must be suitably trained and must have read the instructions given in this manual before proceeding with any operations with it.



After sterilization and before use, check that the first mark on the scale is aligned with the arrow. The instrument must be tested for correct assembly and correct functionality every time it is used.

Important warning

The arm of the torque wrench must not move beyond the end of the scale, as this could lead to inaccurate torque readings.

The torque wrench can also be used as a fixed key, without using the scale, by using the entire handle as a lever. In this case, it must not exceed the torque value of 150 Ncm.

The personnel responsible for the use and maintenance of this instrument must check it for possible signs of deterioration of the tightening, insertion and torque mechanisms.

The single components of the torque wrench are not interchangeable, and it is not possible to use a component from one key to replace a component on another. If any component of the torque wrench is lost, always return the entire instrument to Sweden & Martina S.p.A. for all necessary repairs. Components for the assembly of the torque wrench with control lever are not sold individually. Failure to respect the instructions given may cause aesthetic problems and be damaging for the patient's health.

Cleaning, disinfection, sterilisation and storage of the kit and of the surgical instruments

Attention! All the surgical instruments for dental implants are sold NON-STERILE. Before use, they must be cleaned, disinfected and sterilised according to the following procedure validated by Sweden & Martina. These processes must be performed before use and before each subsequent reuse. Repetition of the processes described in this paragraph has minimal effect on the wear of these devices.

Instruments should always be checked before use to ensure they are in good working order. Any instruments showing signs of wear must be immediately replaced with new devices. It is particularly important to check that the drivers grip properly inside the engagement wells on the heads of the screws to be lifted and tightened with the same. Failure to follow these instructions may cause cross-infection and intraoperative complications.

a. Cleaning

Containers and transport to be used for washing: there are no special requirements. In case of automatic cleaning, use an ultrasound bath with a suitable detergent solution. Use neutral detergents only. Follow the manufacturer's instructions concerning concentrations and washing times. Use demineralised water to prevent the formation of stains and marks. When draining, check the recesses of the devices, holes, etc. to make sure all residues have been completely removed. If necessary, repeat the cycle or clean manually.

When cleaning manually: use a suitable neutral detergent and follow the manufacturer's user instructions. Brush the products with a soft-bristled brush under plenty of running water. Use the brush to apply the detergent to all surfaces. Rinse with distilled water for at least four minutes. Make sure plenty of running water passes through any holes.

For drills with internal irrigation, use the special pins provided with the handpieces to ensure that the irrigation holes are completely clean and free of bone fragments or biological tissues. After rinsing, dry the devices thoroughly and place them inside suitable sterilisation bags. Do not exceed 120 °C when performing a drying cycle in a washing and disinfection appliance.

b. Sterilisation

In a vacuum autoclave, sterilizing as follows:

- autoclave (gravity displacement cycle) at a temperature of 121°C with minimum exposure of 30 minutes and drying cycle of 15 minutes;
- autoclave (dynamic air removal cycle) at the temperature of 132°C with minimum exposure of 4 minutes and drying cycle of 20 minutes.

c. Storage

After sterilisation, the product must remain in the sterilisation bags. The bags should only be opened immediately prior to reuse. In normal conditions, sterilisation bags maintain the sterility of the contents, unless the wrapping is damaged. Therefore, do not use components if the bags in which they were kept are damaged, and resterilise in new bags before using them again. The storage time of products sterilised inside the bags should not exceed that recommended by the manufacturer of the bags. The product must be stored in a cool dry place, away from sunlight, water and sources of heat.

Cleaning, disinfection, sterilisation and storage of the dynamometric ratchet CRI5-KIT

The processes described below must be performed before use and before each subsequent operation. Repetition of the processes described in this paragraph has minimal effect on the wear of the device. The failure to follow these instructions may cause cross infections. Containers and transport to be used for washing: there are no special requirements. As soon as possible after each use, the key must be placed in a container filled with a disinfecting/cleansing solution and covered with a cloth. This prevents the desiccation of the contaminating agents coming from the patient, and dissolves them, thus making cleaning easier and more effective. Completely disassemble the key as shown below:



Completely unscrew the torque adjustment screw and remove the spring inside the handle of the ratchet body. Do not separate the spring from the pin that acts as a stop.



Use the hexagon tip at the bottom of the torque adjustment screw to unscrew and completely remove the connecting screw of the cover from the side marked "OUT". Exert a light pressure in order to avoid damaging the hexagon tip.



After removing the cover, pull out the two components contained inside the ratchet head: the toothed pawl wheel and wheel stop tooth.

In case of manual cleaning, clean the outer and inner surfaces of the instrument mechanically under hot water with a soft bristled brush. Inject hot water using a needleless syringe to wash the hard-to-access holes of the head and the area around the pawl wheel and wheel stop. If necessary, proceed in the same way for the inside of the handle and of the torque adjustment device. Use a suitable neutral detergent and follow the manufacturer's user instructions. Use the brush to apply the detergent to all surfaces. Rinse with distilled water for at least four minutes. Make sure the running water passes abundantly through the passages. In case of automated ultrasound cleaning: use an ultrasound bath with a suitable detergent solution. Use neutral detergents only.

Follow the manufacturer's instructions concerning concentrations and washing times. Use demineralised water to prevent the formation of stains and marks. During this cycle, avoid contact between the pieces because this causes the machined surfaces to deteriorate, and consequently, loss of precision of the torque measurement. When draining, check the recesses of the devices, holes, etc. to make sure all residues have been completely removed. If necessary, repeat the cycle or clean manually.

Please note: blood residues or other deposits reduce the efficacy of the sterilisation process, which is why it is important to clean thoroughly. During cleaning, avoid sprays or jets of liquid and adopt adequate protections. Avoid contact between this instrument and other nickel-plated instruments.

The pieces must be reassembled prior to sterilisation. Dry the parts, lubricate the functional areas lightly and reassemble the key as shown in the figures below.

Too much lubrication may cause the surfaces of the instrument to resurface during sterilisation. Use only the lubricant supplied.



After lubricating the parts shown in the figure, insert the two elements of the ratchet head according to the following sequence: the toothed pawl wheel and then the wheel stop tooth.



Lubricate the contact areas between the tooth of the pawl wheel and the pin of the wheel stop tooth.



Once parts 2 and 3 have been lubricated and inserted in the head of the ratchet body, position the cover and turn the ratchet body from the "OUT" side. Tighten the screw with the hexagon tip of the torque adjustment screw.



Lubricate the spring inside the ratchet handle as shown in the figure. Assemble the torque adjustment screw, making sure the instrument functions properly. Manually activate the pawl wheel.

Sterilisation: in a vacuum autoclave, proceeding as follows:

- Autoclave (Gravity - Displacement Cycles) Temperature of 121°C with a minimum autoclave cycle of 30 minutes and a drying cycle of 15 minutes.

This procedure is important in order to preserve the precision of the instrument within a tolerance of ± 3.5 Ncm. Operate the torque and insertion mechanism to check their proper functioning. Remove any traces of lubricant from the outer surface of the key. Place the device in suitable sterilisation bags. It is recommended to practice the disassembly and reassembly operations, following the instructions.

Cleaning, disinfection, sterilization and storage of the torque wrench with control lever TWL


The processes described below must be performed before use and before each subsequent operation. Repetition of the processes described in this paragraph has minimal effect on the wear of the device. The failure to follow these instructions may cause cross-infections.

a. Cleaning

Containers and transport to be used for washing: there are no special requirements. As soon as possible after each use, the key must be placed in a container filled with a disinfecting/cleansing solution and covered with a cloth.

This prevents the desiccation of the contaminating agents coming from the patient, and dissolves them, thus making cleaning easier and more effective.

Completely disassemble the key as shown below:



Press the driver and remove it from the head of the torque wrench, then remove the head by pressing inside the hollow, and delicately pull it out. The three separate components are now ready for cleaning.

In case of manual cleaning, clean the outer and inner surfaces of the instrument mechanically under hot water with a soft bristled brush. Use a suitable neutral detergent and follow the manufacturer's user instructions. Use the brush to apply the detergent to all surfaces. Rinse with distilled water for at least four minutes. Make sure the running water passes abundantly through the passages. In case of automated ultrasound cleaning: use an ultrasound bath with a suitable detergent solution. Use neutral detergents only. Follow the manufacturer's instructions concerning concentrations and washing times. Use demineralised water to prevent the formation of stains and marks. During this cycle, avoid contact between the pieces because this causes the machined surfaces to deteriorate, and consequently, loss of precision of the torque measurement.

When draining, check the recesses of the devices, holes, etc. to make sure all residues have been completely removed. If necessary, repeat the cycle or clean manually.

Observation: blood residues or other deposits reduce the efficacy of the sterilisation process, which is why it is important to clean thoroughly. During cleaning, avoid sprays or jets of liquid and adopt adequate protections. Avoid contact between this instrument and other nickel-plated instruments.

Components must be reassembled before sterilization.

This procedure is important in order to preserve the precision of the instrument within the following tolerances:

10 Ncm	± 0,75 Ncm
30 Ncm	± 1,5 Ncm
50 Ncm	± 2,5 Ncm
70 Ncm	± 3,5 Ncm
90 Ncm	± 4,5 Ncm



After cleaning, connect the torque wrench head to the body, pushing the components together and rotating them in opposite directions until a click is heard.

Press the driver into the torque wrench until a click is heard. The arrow on the torque wrench head indicates the direction of operation.

Place the device in a suitable sterilization bag. Disassembly and reassembly operations must be carried out following the indications provided here.

b. Sterilization

In a vacuum autoclave, proceeding as follows:

- autoclave (Gravity-Displacement Cycles) at a temperature of 121°C with minimum exposure of 30 minutes and drying cycle of 15 minutes.;
- autoclave (Dynamic –Air – Removal Cycles) Temperature of 132°C with exposition of 4 minutes and a minimum drying cycle of 20 minutes.

c. Storage

After sterilisation, the product must remain in the sterilisation bags. The bags should only be opened immediately prior to reuse. In normal conditions, sterilisation bags maintain the sterility of the contents, unless the wrapping is damaged. Therefore, do not use components if the bags in which they were kept are damaged, and resterilise in new bags before using them again.

The storage time of products sterilised inside the bags should not exceed that recommended by the manufacturer of the bags.

The product must be stored in a cool dry place, away from sunlight, water and sources of heat.

Shorty Drilling Kit

A kit of drills and stops dedicated to Prama cylindrical implants with 6.00 mm height and other short implants (7.00 and 8.50 mm) of Sweden & Martina system is available. Using Shorty drills allows all the available bone to be used for seating the implant since they do not make any overpreparation of the implant site. In addition, the laser markings that report heights from 5.00 mm to 7.00 mm, together with the relative stops, allow safe and rapid preparation.

Compared to traditional drills, the overall length of Shorty drills is shorter (24.85 mm instead of 35.00 mm), making it possible to use them even in case of difficult-to-reach distal sectors or in case of patients with reduced oral opening. In the following pages just the instruments used for the insertion of Prama implants height 6.00 mm are explained.

The kit contains two short parallelism pins, which are very practical in distal sectors



All the instruments contained in the Shorty Drilling Kit are also available as individual spares

The colour codes of the implant diameters facilitate the instrument choice

description	code
Drilling Kit for short implants	ZSHORTY-INT 
Drilling Kit Shorty tray	SHORTY-TRAY-INT 
Kit with 5 spare silicon supports for surgical trays, for drills or instruments with right-angle shank	GROMMET-CA-1 


Important warning

The Shorty Drilling Kit contains only drills and two parallelism pins. It is not a complete surgical kit: to insert implants must be used all the other instruments included in the standard surgical kit (driver, ratchet, screwdrivers, etc.).


Pilot drill

Shorty pilot drill and related stop	
<p>FPS-200</p> 	<p>STOPS4-200-060</p> 

Short guide drills




short guide drills	
<p>FGS-200/300</p> 	<p>FGS-340/425</p> 

Parallelism pin



description	code
<p>parallelism pin for short drills with depth lines at 5.00, 6.00 and 7.00 mm</p>	<p>PPS-2/3</p> 

Final drills and related stops

drills \varnothing	3.80	4.25	5.00
short final drills	FFS-300 	FFS-340 	FFS-425 
stop for preparations H 0.60 mm	STOPS4-300-060 	STOPS4-340-060 	STOPS4-425-060 



- LL:** Total length of the working part, including the tip.
- LS:** Length of the tip. This measurement must be calculated in addition to the length of the preparation hole.

Please note: remember that the drills in the Shorty Drilling Kit do not over-prepare the surgical site. The working lengths include the portion related to the conical tip of the drill.

Drilling Kit Syra Short

The Drilling Kit Syra Short contains the specific drills for the insertion of short conical implants, such as Syra Short and 6.00 mm height Prama RF implants. The drills do not overprepare the implant site and allow to dedicate all the available bone to the implant insertion. Two types of drills are contained in the Drilling Kit Syra Short: one is specific for the preparation in hard bone and the other for the preparation in soft bone, that underprepare the implant site so as to stabilize the implant even in cases of low density bone. The drills have been designed to be used together with the stops, that ensure a more precise and safe preparation. The drills have laser markings at the three different heights of 4.00 mm, 5.00 and 6.00 mm. For the maximum duration of the components, it is recommended to respect the cleaning and sterilisation procedures. In the following pages just the instruments used for the insertion of Prama RF implants height 6.00 mm are explained.

Important warning

The Drilling Kit Syra Short is a set of drills. It is not a complete surgical drill: for the insertion of Prama RF implants the instruments contained in the standard surgical kit are needed (torque wrench, screwdrivers, etc.)



The instruments are made of stainless steel. On the tray the housings of the single components are indicated to simplify the identification of the pieces and their following positioning after cleaning and detersion phases


description	code
Drilling Kit Syra Short for short implants	<p>ZSYRASHORT-INT</p> 
Drilling Kit Syra Short tray	<p>SESHORT-TRAY-INT</p> 

Syra Short pilot drill

The pilot drill SE-FP200 is used to prepare the initial hole for the insertion of the implant. It has a cylindrical morphology and a diameter of 2.00 mm. The drill is recognizable thanks to a white ring and the laser marked code on the shank of the drill. It has a two cutting edges. Besides the laser markings at the heights of 4.00 mm, 5.00 mm and 6.00 mm, the Short pilot drill has also laser markings at standard heights, because it can be also used in distal sectors in case of patients with reduced oral opening.

description	code	
short pilot drill and related stop	SE-FP200 	SE-STOP-FP200-060 stop 6.00 mm for short pilot drill 

LL: Total length of the working part, including the tip



Stop for pilot drill

Stops are devices to be fitted in tip → shank direction on drills suited to receive them. They make it possible to restrict the working length of a drill to a pre-set height. Always check that the stop is inserted at the desired height. Incomplete insertion may reduce the preparation height. Any insertion difficulties can be resolved by loosening the stop tabs slightly, using forceps. It is also recommended to check the retention exerted by the stop, as if retention is too weak the instrument will fall off the drill during operation. In the event of reduced retention capacity, simply tighten the tabs by hand or using forceps.

LL: Total length of the working part, including the tip



Syra Short final drills and related stops

Conical drills are made of stainless steel with high resistance to corrosion and wear. They have a number of cutting edges proportioned to the diameter of the hole, in order to allow a continue and homogeneous cutting movement and a better stability of the instrument during the surgical phases. All this allows to obtain high precision implant preparations, with following easiness in the insertion of the implant.


In order to prepare Prama RF implant site both Syra Short drills for hard bone, which code on the shank is SE-FK*, and Syra Short drills for soft bone, which code on the shank is characterized by the letter “U”, thus the code is SE-FK*-U, must be used.

∅ drills	4.10	5.00	
short final drills	<p>SE-FK410-U drill for preparation in soft bone</p>	<p>SE-FK500 drill for preparation in soft bone</p>	<p>SE-FK500-U drill for preparation in soft bone</p>
stops	<p>SE-STOP-FK410-060</p>	<p>SE-STOP-FK500-060</p>	<p>SE-STOP-FK500-060</p>

Preparation of the implant site

To obtain a three-dimensional view of the bone available, it is recommended to lift a mucoperiosteal flap. As already mentioned previously, pre-operative clinical and radiographic exams play an important role in determining the position and direction according to which the implants will be positioned. In this stage, a surgical stent will be helpful, acting as a guide during the marking of the cortical bone with the Precision drill and in the drilling phase with the 2.00 mm pilot drill.

As a rule a distance of 3.00 mm should be maintained between the perimeter of the implants, and at least 2.00 mm between implants and adjacent natural teeth. The numerous experimental and clinical studies carried out indicate that it is opportune to position the implants more in a lingual or a palatal direction to obtain the best aesthetic results, because this position helps preserve the level of the hard and soft tissues at the crown of the implant. It is also essential to check that the thickness of the residual bone wall at buccal level is not less than 1.00 mm. If the thickness is smaller there is a high risk of bone reabsorption failure and exposure of the spires.



Surgical sequences

The following pages contain information on the drilling sequences for the adequate preparation of all implant types. These procedures come from clinical experience and recommendations taken from numerous studies and clinical protocols for implants of this type. However, it should be remembered that bone types with different densities require different surgical approaches, and the indications below cannot replace the necessary training and knowledge of the doctors, nor their personal experience, which can at times lead to different solutions and indications. The sequences that follow refer to specific bone types. In expansion techniques or in case of regenerative surgery, or when you want to increase the compaction in poor quality bone, the use of drills can be replaced with the relative osteotomes.

Remember to always use drills with stops correctly inserted. Remember that the drills always prepare a hole that is longer than the implant. For the overpreparation dimensions, refer to page 21 for the cylindrical pilot drill, to page 24 for the cylindrical drills to page 28 for the conical drills. The preparations must be non-traumatic and as gradual as possible, and must be executed quickly and precisely. No overheating of the bone should be generated.

It should also be remembered to initially set the surgical micromotor with the correct torque, reduction and rotation values depending on the operation to be performed. In particular:

- the **drill** must be used at the speed indicated in each sequence, with the maximum torque and irrigated copiously with cold sterile physiological solution, better if cooled in a refrigerator;
- the **bone taps** must only be used when indicated in each procedure.

Incorrect insertion of the instruments in the handpiece will cause instrument vibration, eccentric rotation, early wear and shaft buckling. Suitable surgical micromotors only should be used. Micromotors should be checked regularly by their manufacturers, according to the indications given by the same, to prevent potential malfunctions (e.g. axle shifts for transmission shafts, worn or faulty forceps, etc.). Failure to follow the instructions provided may cause surgical problems and damage to the patient's health.

Sweden & Martina distributes **Impla6000** and **Impla7000**, different brushless micromotors for surgical and implant procedures, that combines reliability, high performances and easy to use procedures. Compact, practical and with a basic design, they come with all the requirements for maximum precision and safety.



Impla6000




Impla7000

Surgical sequence for Prama implants

Surgical sequence for Prama Shorty implants H. 6.00 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility.

The graphic sequence is referred to the \varnothing 5.00 mm implant.



	FS-230	FPS-200	FGS-200/300	FFS-300	
LA-ZT-380-060		use up to marking 6.00 mm	use up to half of the last step	use up to marking 6.00 mm	
\varnothing 3.80 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	-	-	-	-
	BONE D4	-	-	-	-
LA-ZT-425-060					
\varnothing 4.25 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	-	-	-	-
	BONE D4	-	-	-	-
LA-ZT-500-060					
\varnothing 5.00 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	-	-	-	-
	BONE D4	-	-	-	-



FFS-340

FGS-340/425

FFS-425

See chart below

EASYC4-EX230-CA

use up to marking 6.00 mm	use up to half of the last step	use up to 6.00 mm	50 Ncm max	50 Ncm max
-	-	-	A-MS-380 (20 rpm)	20 rpm
-	-	-	-	20 rpm
-	-	-	-	-
-	-	-	-	-
			50 Ncm max	50 Ncm max
900 rpm	-	-	A-MS-425 (20 rpm)	20 rpm
900 rpm	-	-	-	20 rpm
-	-	-	-	-
-	-	-	-	-
			50 Ncm max	50 Ncm max
1.100 rpm	1.100 rpm	900 rpm	A-MS-500 (20 rpm)	20 rpm
1.100 rpm	1.100 rpm	900 rpm	-	20 rpm
-	-	-	-	-
-	-	-	-	-


Surgical sequence for Prama implants H. 8.50 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 24 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	LA-ZT-380-085		use up to marking 8.50 mm	use up to marking 8.50 mm	use up to marking 8.50 mm
\varnothing 3.80 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	900 rpm
	BONE D4	900 rpm	900 rpm	osteotome*	osteotome*
LA-ZT-425-085					
\varnothing 4.25 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	900 rpm
	BONE D4	900 rpm	900 rpm	osteotome*	osteotome*
LA-ZT-500-085					
\varnothing 5.00 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	1.100 rpm	900 rpm	900 rpm	900 rpm
	BONE D4	1.100 rpm	900 rpm	osteotome*	osteotome*

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



FFT3-340-LXS

FG-330/425XS

FFT3-425-LXS

See chart below

EASYC4-EX230-CA

use up to marking 8.50 mm	use up to marking 8.50 mm	use up to marking 8.50 mm	50 Ncm max	50 Ncm max
-	-	-	A-MS-380 (20 rpm)	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
900 rpm	-	-	A-MS-425 (20 rpm)	20 rpm
900 rpm	-	-	-	20 rpm
900 rpm	-	-	-	20 rpm
osteotome*	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
1.100 rpm	1.100 rpm	900 rpm	A-MS-500 (20 rpm)	20 rpm
1.100 rpm	1.100 rpm	900 rpm	-	20 rpm
900 rpm	900 rpm	900 rpm	-	20 rpm
osteotome*	osteotome*	osteotome*	-	20 rpm

*All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.


Surgical sequence for Prama implants H. 10.00 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 24 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	LA-ZT-380-100		use up to marking 10.00 mm	use up to marking 10.00 mm	use up to marking 10.00 mm
\varnothing 3.80 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	900 rpm
	BONE D4	900 rpm	900 rpm	osteotome*	osteotome*
\varnothing 4.25 mm	LA-ZT-425-100				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	900 rpm
BONE D4	900 rpm	900 rpm	osteotome*	osteotome*	
\varnothing 5.00 mm	LA-ZT-500-100				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	1.100 rpm	900 rpm	900 rpm	900 rpm
BONE D4	1.100 rpm	900 rpm	osteotome*	osteotome*	

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



FFT3-340-LXS

FG-330/425XS

FFT3-425-LXS

See chart below

EASYC4-EX230-CA

use up to marking 10.00 mm	use up to marking 10.00 mm	use up to marking 10.00 mm	50 Ncm max	50 Ncm max
-	-	-	A-MS-380 (20 rpm)	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
900 rpm	-	-	A-MS-425 (20 rpm)	20 rpm
900 rpm	-	-	-	20 rpm
900 rpm	-	-	-	20 rpm
osteotome*	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
1.100 rpm	1.100 rpm	900 rpm	A-MS-500 (20 rpm)	20 rpm
1.100 rpm	1.100 rpm	900 rpm	-	20 rpm
900 rpm	900 rpm	900 rpm	-	20 rpm
osteotome*	osteotome*	osteotome*	-	20 rpm

*All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.


Surgical sequence for Prama implants H. 11.50 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 24 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	LA-ZT-380-115		use up to marking 11.50 mm	use up to marking 10.00 mm	use up to marking 11.50 mm
\varnothing 3.80 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	900 rpm
	BONE D4	900 rpm	900 rpm	osteotome*	osteotome*
\varnothing 4.25 mm	LA-ZT-425-115				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	900 rpm
BONE D4	900 rpm	900 rpm	osteotome*	osteotome*	
\varnothing 5.00 mm	LA-ZT-500-115				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	1.100 rpm	900 rpm	900 rpm	900 rpm
BONE D4	1.100 rpm	900 rpm	osteotome*	osteotome*	

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



FFT3-340-LXS

FG-330/425XS

FFT3-425-LXS

See chart below

EASYC4-EX230-CA

use up to marking 11.50 mm	use up to marking 10.00 mm	use up to marking 11.50 mm	50 Ncm max	50 Ncm max
-	-	-	A-MS-380 (20 rpm)	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
900 rpm	-	-	A-MS-425 (20 rpm)	20 rpm
900 rpm	-	-	-	20 rpm
900 rpm	-	-	-	20 rpm
osteotome*	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
1.100 rpm	1.100 rpm	900 rpm	A-MS-500 (20 rpm)	20 rpm
1.100 rpm	1.100 rpm	900 rpm	-	20 rpm
900 rpm	900 rpm	900 rpm	-	20 rpm
osteotome*	osteotome*	osteotome*	-	20 rpm

*All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.


Surgical sequence for Prama implants H 13.00 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 24 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	LA-ZT-380-130		use up to marking 13.00 mm	use up to marking 10.00 mm	use up to marking 13.00 mm
\varnothing 3.80 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	900 rpm
	BONE D4	900 rpm	900 rpm	osteotome*	osteotome*
LA-ZT-425-130					
\varnothing 4.25 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	900 rpm
	BONE D4	900 rpm	900 rpm	osteotome*	osteotome*
LA-ZT-500-130					
\varnothing 5.00 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	1.100 rpm	900 rpm	900 rpm	900 rpm
	BONE D4	1.100 rpm	900 rpm	osteotome*	osteotome*

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



FFT3-340-LXS

FG-330/425XS

FFT3-425-LXS

See chart below


EASYC4-EX230-CA

use up to marking 13.00 mm	use up to marking 10.00 mm	use up to marking 13.00 mm	50 Ncm max	50 Ncm max
-	-	-	A-MS-380 (20 rpm)	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
900 rpm	-	-	A-MS-425 (20 rpm)	20 rpm
900 rpm	-	-	-	20 rpm
900 rpm	-	-	-	20 rpm
osteotome*	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
1.100 rpm	1.100 rpm	900 rpm	A-MS-500 (20 rpm)	20 rpm
1.100 rpm	1.100 rpm	900 rpm	-	20 rpm
900 rpm	900 rpm	900 rpm	-	20 rpm
osteotome*	osteotome*	osteotome*	-	20 rpm

*All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.


Surgical sequence for Prama implants H 15.00 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 24 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	FS-230	FPT3-200-LXS	FG-200/280XS	FFT3-300-LXS
\varnothing 3.80 mm	LA-ZT-380-150		use up to marking 15.00 mm	use up to marking 15.00 mm
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	900 rpm	900 rpm	900 rpm
	BONE D4	900 rpm	900 rpm	osteotome*
\varnothing 4.25 mm	LA-ZT-425-150			
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	900 rpm	900 rpm	900 rpm
	BONE D4	900 rpm	900 rpm	osteotome*
\varnothing 5.00 mm	LA-ZT-500-150			
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm
	BONE D3	1.100 rpm	900 rpm	900 rpm
	BONE D4	1.100 rpm	900 rpm	osteotome*

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



FFT3-340-LXS

FG-330/425XS

FFT3-425-LXS

See chart below

EASYC4-EX230-CA


use up to marking 15.00 mm	use up to marking 10.00 mm	use up to marking 15.00 mm	50 Ncm max	50 Ncm max
-	-	-	A-MS-380 (20 rpm)	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
-	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
900 rpm	-	-	A-MS-425 (20 rpm)	20 rpm
900 rpm	-	-	-	20 rpm
900 rpm	-	-	-	20 rpm
osteotome*	-	-	-	20 rpm
			50 Ncm max	50 Ncm max
1.100 rpm	1.100 rpm	900 rpm	A-MS-500 (20 rpm)	20 rpm
1.100 rpm	1.100 rpm	900 rpm	-	20 rpm
900 rpm	900 rpm	900 rpm	-	20 rpm
osteotome*	osteotome*	osteotome*	-	20 rpm

*All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.

Surgical sequence for Prama RF implants

Underpreparation in soft bone

Surgical sequence for Prama RF Shorty implant




ø 4.25 mm	LS-ZT-425-060		use up to marking 6.00 mm	use up to marking 6.00 mm	50 Ncm max
	BONE D2	600 rpm	1.100 rpm	900 rpm	20 rpm
	BONE D3	600 rpm	900 rpm	900 rpm	20 rpm
	BONE D4	600 rpm	900 rpm	900 rpm	20 rpm




ø 5.00 mm	LS-ZT-500-060		use up to marking 6.00 mm	use up to marking 6.00 mm	use up to marking 6.00 mm	50 Ncm max
	BONE D2	600 rpm	1.100 rpm	900 rpm	900 rpm	20 rpm
	BONE D3	600 rpm	900 rpm	900 rpm	900 rpm	20 rpm
	BONE D4	600 rpm	900 rpm	900 rpm	900 rpm	20 rpm

Preparation in hard bone

Surgical sequence for Prama RF Shorty implant




ø 4.25 mm	LS-ZT-425-060		use up to marking 6.00 mm	use up to marking 6.00 mm	use up to marking 6.00 mm	50 Ncm max
	BONE D1	800 rpm	1.100 rpm	900 rpm	900 rpm	20 rpm



ø 5.00 mm	LS-ZT-500-060		use up to marking 6.00 mm	50 Ncm max			
	BONE D1	800 rpm	1.100 rpm	900 rpm	900 rpm	900 rpm	20 rpm

Surgical sequence for Prama RF implants H. 8.50 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 28 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	FS-230	FPT3-200-LXS	FG-200/280XS	SH-FK380-085	
\varnothing 3.80 mm	LS-ZT-380-085		use up to marking 8.50 mm	use up to marking 8.50 mm	
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-
\varnothing 4.25 mm	LS-ZT-425-085				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-
\varnothing 5.00 mm	LS-ZT-500-085				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



SH-FK425-085

SH-FK500-085

See chart below

EASYC4-EX230-CA

		50 Ncm max	50 Ncm max
-	-	SH-MS-380-CA (20 rpm)	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm

		50 Ncm max	50 Ncm max
900 rpm	-	SH-MS-425-CA (20 rpm)	20 rpm
900 rpm	-	-	20 rpm
800 rpm	-	-	20 rpm
-	-	-	20 rpm

		50 Ncm max	50 Ncm max
900 rpm	900 rpm	SH-MS-500-CA (20 rpm)	20 rpm
900 rpm	900 rpm	-	20 rpm
800 rpm	800 rpm	-	20 rpm
-	-	-	20 rpm

*All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.

Surgical sequence for Prama RF implants H. 10.00 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 28 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	LS-ZT-380-100		use up to marking 10.00 mm	use up to marking 10.00 mm	
\varnothing 3.80 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-
\varnothing 4.25 mm	LS-ZT-425-100				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
BONE D4	900 rpm	osteotome*	-	-	
\varnothing 5.00 mm	LS-ZT-500-100				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
BONE D4	900 rpm	osteotome*	-	-	

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



SH-FK425-100

SH-FK500-100

See chart below

EASYC4-EX230-CA

		50 Ncm max	50 Ncm max
-	-	SH-MS-380-CA (20 rpm)	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm


		50 Ncm max	50 Ncm max
900 rpm	-	SH-MS-425-CA (20 rpm)	20 rpm
900 rpm	-	-	20 rpm
800 rpm	-	-	20 rpm
-	-	-	20 rpm

		50 Ncm max	50 Ncm max
900 rpm	900 rpm	SH-MS-500-CA (20 rpm)	20 rpm
900 rpm	900 rpm	-	20 rpm
800 rpm	800 rpm	-	20 rpm
-	-	-	20 rpm

**All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.*

Surgical sequence for Prama RF implants H. 11.50 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 28 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	LS-ZT-380-115		use up to marking 11.50 mm	use up to marking 10.00 mm	
\varnothing 3.80 mm	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-
\varnothing 4.25 mm	LS-ZT-425-115				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
\varnothing 5.00 mm	LS-ZT-500-115				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



SH-FK425-115

SH-FK500-115

See chart below

EASYC4-EX230-CA

		50 Ncm max	50 Ncm max
-	-	SH-MS-380-CA (20 rpm)	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm


		50 Ncm max	50 Ncm max
900 rpm	-	SH-MS-425-CA (20 rpm)	20 rpm
900 rpm	-	-	20 rpm
800 rpm	-	-	20 rpm
-	-	-	20 rpm

		50 Ncm max	50 Ncm max
900 rpm	900 rpm	SH-MS-500-CA (20 rpm)	20 rpm
900 rpm	900 rpm	-	20 rpm
800 rpm	800 rpm	-	20 rpm
-	-	-	20 rpm

*All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.

Surgical sequence for Prama RF implants H. 13.00 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 28 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	FS-230	FPT3-200-LXS	FG-200/280XS	SH-FK380-130	
\varnothing 3.80 mm	LS-ZT-380-130		use up to marking 13.00 mm	use up to marking 10.00 mm	
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-
\varnothing 4.25 mm	LS-ZT-425-130				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-
\varnothing 5.00 mm	LS-ZT-500-130				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



SH-FK425-130

SH-FK500-130

See chart below

EASYC4-EX230-CA

		50 Ncm max	50 Ncm max
-	-	SH-MS-380-CA (20 rpm)	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm

		50 Ncm max	50 Ncm max
900 rpm	-	SH-MS-425-CA (20 rpm)	20 rpm
900 rpm	-	-	20 rpm
800 rpm	-	-	20 rpm
-	-	-	20 rpm

		50 Ncm max	50 Ncm max
900 rpm	900 rpm	SH-MS-500-CA (20 rpm)	20 rpm
900 rpm	900 rpm	-	20 rpm
800 rpm	800 rpm	-	20 rpm
-	-	-	20 rpm

*All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.

Surgical sequence for Prama RF implants H. 15.00 mm

The use of the stops depends on the clinician. Their use is recommended, most of all in the case of poor intraoperative visibility. Please remember that the drills over-prepare the length for a measurement reported in the table on page 21 (for the pilot drill) and 28 (for the final drills). The graphic sequence is referred to the \varnothing 5.00 mm implant.



	FS-230	FPT3-200-LXS	FG-200/280XS	SH-FK380-150	
\varnothing 3.80 mm	LS-ZT-380-150		use up to marking 15.00 mm	use up to marking 10.00 mm	
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-
\varnothing 4.25 mm	LS-ZT-425-150				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-
\varnothing 5.00 mm	LS-ZT-500-150				
	BONE D1	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D2	1.100 rpm	1.100 rpm	1.100 rpm	900 rpm
	BONE D3	900 rpm	900 rpm	900 rpm	800 rpm
	BONE D4	900 rpm	osteotome*	-	-

In the case of surgeries in the distal sectors or in case of poor oral opening of the patient, drills of reduced length are available, to be used without stop. For further details see page 34.



SH-FK425-150

SH-FK500-150

See chart below

EASYC4-EX230-CA

		50 Ncm max	50 Ncm max
-	-	SH-MS-380-CA (20 rpm)	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm
-	-	-	20 rpm

		50 Ncm max	50 Ncm max
900 rpm	-	SH-MS-425-CA (20 rpm)	20 rpm
900 rpm	-	-	20 rpm
800 rpm	-	-	20 rpm
-	-	-	20 rpm

		50 Ncm max	50 Ncm max
900 rpm	900 rpm	SH-MS-500-CA (20 rpm)	20 rpm
900 rpm	900 rpm	-	20 rpm
800 rpm	800 rpm	-	20 rpm
-	-	-	20 rpm

**All the osteotomes have to be used at the reference notch of the implant to insert.
For further details see page 36.*

Implant insertion

1 Use the patient label found inside the pack for the patient's medical file and apply it on the Dental Card: this will make it easier to record the patient's treatment plan and will keep a trace of the batch used.



2 Open the blister and place the vial contained in it on a sterile surface (i.e. on a disposable towel or sterile cloth) next to the operating field.




3 Immediately before inserting it into the oral cavity, remove the blue cap of the vial, making sure not to remove the transparent cap containing the surgical cover screw. The implant holding cylinder inside the vial and the surgical cover screw are coloured according to a colour code that allows the rapid identification of the implant diameter.



Standard procedure

When the vial is opened the mounter is presented with the hexagon ready to be engaged. The implant may be picked up using the dedicated driver and then screwed mechanically in place with the aid of a suitable surgical micromotor with torque control set at a screwing speed of 20 rpm and max torque 70 Ncm. At the moment this value is the maximum that can be reached by the micromotors on the market. The driver has been tested up to 70 Ncm and has not presented any deformations or failures. Instruments with torque control, both mechanical and normal, are regularly calibrated with a suitable calibrated instrument.



Phase after inserting the implant

Healing times

It is essential to respect the healing times recommended in implant surgery and to check periodically the state of evolution of osseointegration, with x-rays, the stadium of the osseointegration. The preliminary healing times of an implant are influenced by the quality of the receiving bone. In the case of immediate load, consider the warnings reported on pages 6-7. In the case of a deferred loading, in order to minimize the discomfort conditioned by the biological times for the osseointegration, the use of mobile temporary prosthesis has to be carried out with prudence, widely unloading the prosthesis.

After the healing the screws are removed from the implants. If the right angle driver is used, the surgical micromotor must be set with the following parameters: 20 rpm and torque 10 Ncm. After that, depending on the protocol adopted, proceed with the adaptation of the profiles of the tissues with a proper temporary or with proper healing abutments.

It is recommended to secure the healing abutments manually or at any rate with a torque no greater than 10 Ncm.

Intraoperative removal of implants if necessary

If a previously inserted implant needs to be removed, this can be done by directly engaging the hexagonal driver connection of the implant. Accurately clean away blood and any other residues produced during insertion from the implant socket by irrigating the site.



Take the BC-EX230 or BL-EX230 driver not included in the surgical kit and insert the hexagonal tip of the driver inside the implant connection, being very careful that the instrument is on-axis with the implant and that it completely and closely engages the internal hexagon.



Block the head of the CRI5-KIT ratchet or of the TWL key and connect it with the hexagonal tip of the driver making sure that the laser-etched arrow on the ratchet head indicates an anticlockwise direction, and move it in this direction while keeping the driver/ratchet assembly on-axis with the index finger. It is recommended to apply a higher torque than the one applied during the insertion phase. Once it has been unscrewed pick up the removed implant using sterile forceps.



Maintenance of the prosthesis

Some implant restoration-related complications are reported in the literature. These complications may lead to a loss of osseointegration and implant failure. Correct maintenance by the patient, good home dental care and regular sessions with a professional hygienist increase the device's service life. Complications such as the pull-out of screws that fasten the restoration to the implants or bone reabsorption causing the loss of the mucosal resting surface in patients with removable restorations can be easily prevented with regular check-ups. If post or prosthetic connecting screws are needed, these operations must be performed by the practitioner using suitable devices with torque tightening control. The calibration of these devices should be checked regularly. In the event of complications of this kind, patients should contact their practitioner as soon as possible, so that the restoration can be repaired and functionality restored.

A delay in contacting the doctor may lead to the fracture of the connecting screw or of the prosthesis, in the first case, and to implant failure in the second case, which could impair the rehabilitative result. Practitioners must make this clear to their patients. Complications can be of a biological nature (loss of integration) or mechanical nature (fracture of a component due to overloading). If there are no complications, duration depends on the devices and the whole restoration system depends on mechanical resistance in relation to the fatigue accumulated by the device.

Responsibility for defective products and warranty terms

Optimal patient care and attention to their needs are necessary conditions for the success of implantation procedures and, therefore, patients must be carefully selected and informed of the associated risks and obligations connected with the treatment and encouraged to cooperate with the odontologist in the interests of the success of the same treatment. The patient must, therefore, maintain good hygiene, which should be confirmed during check-up appointments, guaranteed and recorded and the practitioners instructions and orders shall be observed. Sweden & Martina offers unlimited lifetime warranty for defects as long as the faulty piece is identified by the article code and batch number and returned within the validity period of the warranty. The warranty terms are available on the website www.sweden-martina.com.

Disposal

If removed from the oral cavity due to biological or mechanical failure, the implant fixtures must be disposed of as biological waste. The surgical instruments are made of small components, mostly metal. They may be disposed of as such. If dirty, they must be disposed of as biological waste. In general, the local regulations apply.

Material composition

The materials used for manufacturing the surgical instruments illustrated in this manual were selected based on the properties indicated for their intended use according to directive 93/42, implemented in Italy with Law 46/97, Annex I – Essential Requirements, point 7.1.

Implants

The implants are made of Gr. 4 commercially pure titanium, in compliance with harmonized standards. Although very rare, titanium allergy is possible. Patients should therefore always be asked whether they have allergies of this type.

The characteristics of the Gr. 4 titanium used are listed below.

Grade 4 titanium (Cold worked)*

chemical composition	maximum allowed values (%)	tolerance
nitrogen	0.05	+/- 0.02
carbon	0.10	+/- 0.02
hydrogen	0.015	+/- 0.002
iron	0.25	+/- 0.10 (%<0.25) +/- 0.15 (%>0.25)
oxygen	0.20	+/- 0.02 (%<0.20) +/- 0.03 (%>0.20)
titanium	remainder	-

* This technical information complies with the express specifications of the regulations in force on the use of Gr. 4 titanium in implantology:

- ASTM F67 – 13(2017): Standard Specification for unalloyed titanium, for surgical implant applications.
- ISO 5832-2:2018: Implant for surgery - Metallic materials - Part 2: Unalloyed titanium.

Important note: the use of cold worked Gr. 4 titanium bars for the production of Sweden & Martina implants allows the exploitation of mechanical characteristics higher than those required by applicable standards. Furthermore, the excellent results documented during 20 years of clinical experience corroborate the choice of the coldworking production process and of ZirTi surface treatments, which express and enhance the raw material potential selected by Sweden & Martina.

Surgical instruments

Depending on the type of component, surgical instruments are made of:

- Gr. 5 titanium
- 1.4197 steel
- 1.4542 steel
- 1.4305 (AISI 630) steel
- 1.4108 (AISI 303) steel
- 1.4108 steel
- 1.4112 steel

Patients must be asked if they are allergic to any of the materials used.

Identification of the manufacturer

The manufacturer of Premium One implants and of the respective surgical instruments is:

Sweden & Martina

Via Veneto 10 - 35020 Due Carrare (Padua) - Italy

Tel. +39 049.9124300 - Fax + 39 049.9124290

e-mail: info@sweden-martina.com

www.sweden-martina.com

In accordance with Directive 93/42/EEC implemented in Italy with L.D. 46/97 of 26/03/97, Annex IX, Sweden & Martina identifies the risk class of these products as indicated in Table 01. Even though they can be used with all patients who have suitable therapeutic indications, dental implants and the respective surgical instruments must be used only by professional dentists or surgeons with the necessary qualifications and training.

Table 01 - Risk classes

device	directive 93/42	packing	rule Annex IX	risk class
Implant fixtures for dental use, belonging to the Prama implant system	Implantable devices intended for long-term use (over 30 days)	Single-use and sterile package, fixture complete with surgical cover screw	8	IIb
Surgical cover screws	Implantable devices intended for long-term use (over 30 days)	Sold in packages complete with the respective fixtures or sold individually (single-use and sterile packages)	8	IIb
Complete surgical kits	Reusable surgical instruments	Sold in NON sterile packages	6	IIa
Radel instrument trays and x-ray templates	Non invasive medical devices	Sold in NON sterile packages	1	I
Surgical drills (precision, conical, cylindrical, for distal use, countersinks, bone profilers) and Drill extensions, Drill stops, Bone taps, Drivers and Drivers/Screwdrivers intended to be used with a micromotor	Reusable invasive surgical instruments for temporary use (for less than 60 minutes at a time)	Sold in NON sterile packages	6	IIa
Osteotomes/Bone Expanders, Drivers/ Screwdrivers, Bone taps, Drivers, Hex drivers, Hand knobs, Depth gauges, Parallelism pins and Stents	Reusable surgical instruments for temporary use (for less than 60 minutes at a time), not intended to be connected to an active medical device	Sold in NON sterile packages	6	I

Key to symbols used on implant packs:

description	symbol
Caution! See instruction for use	
Batch number	
Code	
Manufacturer	
Consult instruction for use	
CE conformity mark for class IIa and IIb products	
American federal law restricts this device to sale by or by order of a professional practitioner	
Do not resterilize	
Disposable product, do not reuse	
Do not use if the packaging is damaged	
Sterilized with ionizing radiation	
Expiry date after which the product must not be used	

Key to symbols used on surgical instrument packs:

description	symbol
Caution! See instruction for use	
Batch number	
Code	
Manufacturer	
Consult instruction for use	
CE conformity mark for class IIa and IIb products	
CE conformity mark for class I products	
American federal law restricts this device to sale by or by order of a professional practitioner	
Non-sterile product	

Key to symbols used on prosthesis packs:

description	symbol
Caution! See instruction for use	
Batch number	
Code	
Manufacturer	
Consult instruction for use	
CE conformity mark for class IIa and IIb products	
CE conformity mark for class I products	
American federal law restricts this device to sale by or by order of a professional practitioner	
Disposable product, do not reuse	
Non-sterile product	

THIS MANUAL WAS LAST UPDATED IN JULY 2018.

The medical devices addressed by this manual have been designed and manufactured in accordance with the most recent directives and harmonized standards applicable to the materials used, production processes, the information supplied and packaging.

Bibliography on Sweden & Martina implants since 2013

For older publications please consult
"Scientifica - Implantology Bibliographic Review"

- Gandolfi M.G., Siboni F., Piattelli A., Prati C.; Nano-topography, microchemical properties and calcium phosphates nucleation of premium implant surfaces; 30th Annual Congress, American Academy of Osseointegration, San Francisco, 12-14 March 2015, Poster Id 2088727
- Caneva M., Lang N.P., Calvo Guirado J.L., Spriano A.M., Iezzi G., Botticelli D.; Bone healing at bicortically installed implants with different surface configurations. an experimental study in rabbits; *Clinical Oral Implant Research*, 2015; 26:293-299
doi: 10.1111/clr.12475
- Beolchini M, Lang N.P., Gomez Moreno G., Iezzi G., Botticelli D., Calvo Guirado J.L.; Bone healing at implants with different surface configurations: an experimental study in dogs, *Clinical Oral Implant Research*, 2015; 00:1-7, doi: 10.1111/clr.12562
- Baffone G., Lang N.P., Pantani F., Favero G., Ferri M., Botticelli D.; Hard and soft tissue changes around implants installed in regular-sized and reduced alveolar bony ridges. An experimental study in dogs. *Clinical Oral Implant Research*, 2015; 26:96-101; doi: 10.1111/clr.12306
- Bengazi F., Lang N.P., Caroprese M., Velez J.U., Favero V., Botticelli D.; Dimensional changes in soft tissues around dental implants following free gingival grafting: an experimental study in dogs; *Clinical Oral Implant Research*, 2015; 26:176-82; doi: 10.1111/clr.12280
- Morelli F., Lang N.P., Bengazi F., Baffone D., Vila Morales C.D., Botticelli D.; Influence of bone marrow on osseointegration in long bones: an experimental study in sheep; *Clinical Oral Implant Research*, 2015; 26:300-306; doi: 10.1111/clr.12487
- Mainetti T., Lang N., Bengazi F., Sbricoli L., Soto Cantero L., Botticelli D.; Immediate loading of implants installed in a healed alveolar bony ridge or immediately after tooth extraction: an experimental study in dogs; *Clinical Oral Implant Research*, 2015; 26:435-441; doi: 10.1111/clr.12389
- Borgia V., Alfonsi F., Toti P., Tonelli P., Covani U., Barone A.; Immediate restoration of post-extraction implants. a 7 years prospective single cohort study; 30th Annual Congress, American Academy of Osseointegration, San Francisco, 12-14 March 2015, Poster
- Guazzotti P.P.; Carico immediato di impianti post estrattivi: presentazione di un caso clinico full-arch; *Doctor OS*, 2015; XXVI, 01
- Agustín Panadero R., Serra Pastor B., Chust López C., Fons Font A., Ferreira A.; Immediate placement of single implant simultaneously with immediate loading in a fresh socket associated to periapical infection: a clinical case report; *Journal of Clinical and Experimental Dentistry* 2015;7(1):e175-9
- Crespi R., Bruschi G. B., Gastaldi G., Cappare P., Gherlone E.F.; Immediate loaded implants in split-crest procedure; *Clinical Implant Dentistry and Related Research*, Article first published online: 17 MAR 2015; DOI: 10.1111/cid.12316
- Peñarrocha Oltra D., Covani U., Peñarrocha Diago M., Peñarrocha Diago M.A.; Immediate versus conventional loading with fixed full-arch prostheses in mandibles with failing dentition: a prospective controlled study; *The International Journal of Oral & Maxillofacial Implants* 2015;30:427-434; doi: 10.11607/jomi.3534
- Canullo L., Peñarrocha Oltra D., Covani U., Botticelli D., Serino G., Peñarrocha Diago M.; Clinical and microbiological findings in patients with peri-implantitis: a cross-sectional study; *Clinical Oral Implants Research* 2015; 00:1-7; doi: 10.1111/clr.12557
- Requena Gómez E., Cervantes Haro M.N., Aragonese Lamas J.M.; ¿Es la cirugía guiada junto a la carga inmediata una técnica predecible? a propósito de un caso clínico; *Numeri Uno* 2015; 04: 16-19
- Peñarrocha Oltra D., Covani U., Peñarrocha Diago M., Peñarrocha Diago M.A.; Immediate versus conventional loading for the maxilla with implants placed into fresh and healed extraction sites to support a full-arch fixed prosthesis: nonrandomized controlled clinical study; *The International Journal of Oral & Maxillofacial Implants* 2015;30:427-434; doi: 10.11607/jomi.3534
- Bruschi G.B., Crespi R., Cappare P., Grande N., Bruschi E., Gherlone E.; Radiographic evaluation of crestal bone levels of delayed implants at medium-term follow-up; *The International Journal of Oral & Maxillofacial Implants* 2014;29:441-447
doi: 10.11607/jomi.3254
- Prati C., Zamparini F., Ciulla A., Buonavoglia A., Gatto M.R., Piattelli A., Gandolfi M.G.; Evaluation of marginal bone level of premium implants; IADR General Session, Boston 11-14 Marzo 2015, Poster
- Canullo L., Peñarrocha Oltra D., Soldini C., Mazzocco F., Peñarrocha Diago M., Covani U.; Microbiological assessment of the implant-abutment interface in different connections: cross-sectional study after 5 years of functional loading; *Clinical Oral Implantology*, 2015; 26:426-434, doi: 10.1111/clr.12383
- Kern J.S., Kern T., Wolfart S., Heussen N.; Review - a systematic review and meta-analysis of removable and fixed implant-supported prostheses in edentulous jaws: post-loading implant loss; *Clinical Oral Implants Research* 2015; 00:1-22 ; doi: 10.1111/clr.12531
- Martín Anciburo M.A.; Rehabilitación unitaria implantosoportada utilizando la técnica B.O.P.T. ,*Numeri Uno* 2015; 04:11-14
- Agustín Panadero R., Serra Pastor B., Roig Vanaclocha A., Román Rodríguez J.L., Fons Font A.; Mechanical behavior of provisional implant prosthetic abutments; *Medicina Oral, Patología Oral y Cirugía Bucal* 2015; 20(1):e94-102
- Crespi R., Cappare P., Polizzi E.M., Gherlone E.F. ; Tissue remodeling after bone expansion in grafted and ungrafted sockets *The International Journal of Oral & Maxillofacial Implants*, 2014;29:699-704; doi: 10.11607/jomi.3535
- Negri B., López Marí M., Maté Sánchez de Val J.E., Iezzi G., Bravo González L.A., Calvo Guirado J.L.; Biological width formation to immediate implants placed at different level in relation to the crestal bone: an experimental study in dogs; *Clinical Oral Implant Research*, 2014; 00:1-11 ;doi: 10.1111/clr.12345
- Esposito M., Ardebili Y., Worthington H.V.; Interventions for replacing missing teeth: different types of dental implants (review); *Cochrane database of systematic reviews*, 2014;22:7; doi: 10.1002/14651858.CD003815.pub4.
- Canullo L., Peñarrocha Oltra D., Peñarrocha Diago M., Rocio A.G., Peñarrocha Diago M.A.; Piezoelectric vs. conventional drilling in implant site preparation: pilot controlled randomized clinical trial with crossover design; *Clinical Oral Implants Research* 2014; 25:1336-43; doi: 10.1111/clr.12278
- Lumetti S., Di Blasio A., Manfredi E., Ghiacci G., Toffoli A., Bonanini M., Macaluso G.M., Galli C.; Implant surface microtopography affects cell the pattern of cell growth, cell-to-cell contacts and the expression of connexin 43; *Clinical Oral Implant Research*, 2014; 25 Suppl 10:222
- Negri M., Galli C., Smerieri A., Macaluso G.M., Manfredi E., Ghiacci G., Toffoli A., Bonanini M., Lumetti S.; The effect of age, gender and insertion site on marginal bone loss around endosseous implants: results from a 3-year trial with premium implant system; *BioMed research International*, 2014; Article ID 369051; 7; doi: org/10.1155/2014/369051
- Quaranta A., Andreana S., Pompa G., Procaccini M.; Active implant peri-apical lesion: a case report treated via guided bone regeneration with a 5-year clinical and radiographic follow-up ; *Journal of Oral Implantology* 2014;40:313-319; doi: 10.1563/AAID-JOI-D-11-00214
- Bowen Antolín A., Ariño B., Arlandi Garrido M.; Regeneración ósea perimplantaria con fosfato de calcio bifásico y ácido poliláctico; *Gaceta Dental*, 2014, 260(7): 174-186
- Mainetti T., Lang N.P., Bengazi F., Favero V., Soto Cantero L., Botticelli D.; Sequential healing at implants installed immediately into extraction sockets. An experimental study in dogs; *Clinical Oral Implant Research*, 2014; 00:1-9; doi: 10.1111/clr.12533
- Covani U., Marconcini S., Ferrini F., Gelpi F., Finotti M., Barone A.; Post-traumatic use of dental implants immediately after tooth extraction - clinical study, *The Journal of Craniofacial Surgery*, 2014; 25:796-798; doi 10.1097/SCS.0000000000000522
- Engelhardt S., Papacosta S., Rathe F., Ozen J., Jansen J.A., Junker R.; Annual failure rates and marginal bone-level changes of immediate compared to conventional loading of dental implants. a systematic review of the literature and meta-analysis; *Clinical Oral Implants Research* 2014;00:1-17; doi: 10.1111/clr.12363
- Romanos G.R., Javed F.; Platform switching minimises crestal bone loss around dental implants: truth or myth? *Journal of Oral Rehabilitation*, 2014; 41:700-708; doi: 10.1111/joor.12189
- Strietzel F.P., Neumann K., Hertel M.; Review article: impact of platform switching on marginal peri-implant bone-level changes. a systematic review and meta-analysis. *Clinical Oral Implant Research*, 2014; 00:1-16; doi: 10.1111/clr.12339
- Kinaia B.M., Shah M., Neely A.L., Goodies H.E.; Crestal bone level changes around immediately placed implants: a systematic review and meta-analyses with at least 12 months' follow-up after functional loading; *Journal of Periodontology*, 2014; 85:1537-48; doi: 10.1902/jop.2014.130722. Epub 2014 May 2

- Covani U., Canullo L., Toti P., Alfonsi F., Barone A.; Tissue stability of implants placed in fresh extraction sockets: a 5-year prospective single-cohort study ; *Journal of Periodontology*, 2014; 85:e323-332; doi: 10.1902/jop.2014.140175. Epub 2014 May 16.
- D'Ercole S., Tripodi D., Marzo G., Bernardi S., Continenza M.A., Piatelli A., Iaculli F., Mummolo S.; Microleakage of bacteria in different implant-abutment assemblies: an in vitro study ; *Journal of Applied Biomaterial and Functional Materials*, 2014, accepted June 12; doi: 105301/jabfm.5000214
- Peñarrocha Oltra D., Rossetti P.H., Covani U., Galluccio F., Canullo L.; Microbial leakage at the implant/abutment connection due to implant insertion maneuvers: cross-sectional study 5 years post loading in healthy patients; *Journal of Oral Implantology*, 2014; 23 [Epub ahead of print]
- Maiorana C., Farronato D., Pieroni S., Ciccì M., Andreoni D., Santoro F.; A four-year survival rate multicenter prospective clinical study on 377 implants: correlations between implant insertion torque, diameter and bone quality; *Journal of Oral Implantology* 2014;11 [Epub ahead of print]
- Crespi R., Bruschi G.B., Capparè P., Gherlone E.; The utility of the electric mallet; *The Journal of Craniofacial Surgery*, 2014;25;793-795; doi 10.1097/SCS.0000000000000523
- Schirripa G., Schirripa F.; Carico immediato; *Numeri Uno*, 2014, 19, 22-24
- Csonka M.; Trattamento implantologico delle creste sottili: split crest o gbr? ; *Numeri Uno*, 19: 12-14, 2014
- Machín Muñiz A.; Regeneración ósea y gingival en implantes inmediatos post-extracción; *Numeri Uno* 2014; 01: 20-21
- Peñarrocha Oltra D., Peñarrocha Diago M.A., Canullo L., Covani U., Peñarrocha Diago M.; Patient-reported outcomes of immediate versus conventional loading with fixed full-arch prostheses in the maxilla: a nonrandomized controlled prospective study; *The International Journal of Oral & Maxillofacial Implants*, 2014;29:690-698; doi: 10.11607/jomi.3516
- Baldi D., Colombo J., Pera P., Hauschild U.; Una tecnica minimamente invasiva: implantologia con utilizzo di impianti a diametro ridotto e tecniche cad cam per una provvisorizzazione a lungo termine; *Numeri Uno*, 2014;18: 6-9
- Calesini G., Zarone F., Sorrentino R., Micarelli C., Fabianelli A., Papacchini F., Gherlone E.; Effect of 2 impression techniques on the dimensional accuracy of working implant prosthesis models: an in vitro study; *Journal of Craniofacial Surgery* 2014;25:822-827
- Pellicer Chover H., Peñarrocha Oltra D., Bagán L., Fichy Fernandez A.J., Canullo L., Peñarrocha Diago M.; Single-blind randomized clinical trial to evaluate clinical and radiological outcomes after one year of immediate versus delayed implant placement supporting full-arch prostheses; *Medicina Oral Patología Oral y Cirugía Bucal*, 2014; 19: e295-301
- Morandini E.; La precisione nel cr.co. laser sinterizzato rivestito in ceramica parte 2; *NumeriUno*, 2014;18: 16-19
- De Santis E., Lang N.P., Favero G., Beolchini M., Morelli F., Botticelli D.; Healing at mandibular block-grafted sites. an experimental study in dogs; *Clinical Oral Implant Research*, 2014; 00:1-7; doi: 10.1111/clr.12434
- Cocchetto R.; Improved cementation technique for implant restorations to avoid peri-implant cement remnants: clinical and microscopical evaluation with two different abutment design; *Clinical Oral Implants Research* 2014; 25(Suppl. 10); Doi 10.1111 clr.12458_94
- J. Viña Almunia; Microbial colonization of the implant connection with cemented versus screw-retained suprastructures *Clinical Oral Implants Research*, 2014; 25; DOI 10.1111/clr.12458_91
- Ciccì M., Bramanti E., Maticena G., Guglielmino E., Risitano G.; Fem evaluation of cemented-retained versus screw-retained dental implant single-tooth crown prosthesis ; *International Journal of Clinical and Experimental Medicine* 2014; 7(4):817-825
- Vischia F., Roncoroni F.; Ortodonzia protesica mediante tecnica B.O.P.T. ; *Numeri Uno*, 2014;19:19-21
- Loi I.; Tecnica B.O.P.T. su denti e impianti per la riabilitazione di un'arcata completa; *Numeri Uno*, 2014;18:21-22
- Vedove F.; Riabilitazione di elemento singolo in zona estetica con impianto Prama; *Numeri Uno*, 2014;20:18-19
- Gorni F.; Riabilitazione di elemento singolo in zona estetica con impianto Prama RF; *Numeri Uno*, 2014;20:16-17
- Andreoni D.; Riabilitazione di elemento singolo in posizione 4.6 con impianto Prama; *Numeri Uno*, 2014; 20: 20-21
- Sandri L.P.; Utilizzo clinico dei nuovi impianti Prama: inserimento e riabilitazione con un singolo impianto; *Numeri Uno* 2014; 20:22-24
- Loi I.; Riabilitazione impianto-protesica di elemento incisivo frontale con impianto Prama; *Numeri Uno*, 2014; 20:12-13
- Loi I.; Riabilitazione implantoprotesica di ponte distale con impianti Prama ;*Numeri Uno*, 2014; 20:14-15
- Canullo L., Peñarrocha Oltra D., Marchionni S., Bagán L., Peñarrocha Diago M.A., Micarelli C.; Soft tissue cell adhesion to titanium abutments after different cleaning procedures: preliminary results of a randomized clinical trial; *Medicina Oral, Patología Oral y Cirugía Bucal* 2014;19(2):e177-83
- Canullo L., Micarelli C., Bettazzoni L., Magnelli A., Baldissara P.; Shear bond strength of veneering porcelain to zirconia after argon plasma treatment; *The International Journal of Prosthodontics* 2014;27(2):137-139; doi: 10.11607/ijp.3722
- Canullo L., Micarelli C., Bettazzoni L., Koçi B., Baldissara P.; Zirconia-composite bonding after plasma of argon treatment; *The International Journal of Prosthodontics* 2014; 27:267-269; doi: 10.11607/ijp.3686
- Marchetti E., Ratta S., Mummolo S., Tecco S., Pecci R., Bedini R., Marzo G.; Evaluation of an endosseous oral implant system according to uni en iso 14801 fatigue test protocol; *Implant Dentistry*, 2014, Early View in ahead of print; doi: 10.1097/id.151
- Crespi R., Capparè P., Gastaldi G., Gherlone E.F.; Immediate occlusal loading of full-arch rehabilitations: screw-retained versus cement-retained prosthesis. an 8 year clinical evaluation; *International Journal of Oral & Maxillofacial Implants* 2014;29:1406-1411; doi: 10.11607/jomi.3746
- Peñarrocha Oltra D., Candel Martí M.E., Peñarrocha Diago M., Agustín-Panadero R., Canullo L., Peñarrocha Diago M.A.; The horizontal denture: a prosthodontic alternative for patients with severe maxillary atrophy. a technical note; *Journal of Oral Implantology* 2014; 8 [Epub ahead of print]
- Gaspari L.; Tecnica conometrica con provvisorio elettrosaldato per carico immediato; *Italian Dental Journal* 2014; 29, agosto
- Gaspari L.; Implantoprotesi conometrica elettrosaldato chairside a carico immediato - caso clinico; *Numeri Uno* 2014;18:12-14
- Pradies Ramiro G., Abad Coronel C., García Martínez I., Ferreiroa Navarro A.; Impresiones fiables: dos propuestas para un mismo objetivo; *Numeri Uno* 2014; 01:6-9
- Beolchini M., Lang N.L., Ricci E., Bengazi F., Garcia Triana B., Botticelli D.; Influence on alveolar resorption of the buccal bony plate width in the edentulous ridge expansion (e.r.e.) – an experimental study in the dog; *Clinical Oral Implant Research*, 2013; 00:1-6;doi: 10.1111/clr.12308
- Petrillo N.; Carico immediato full-arch mascellare e mandibolare: un nuovo approccio chirurgico e protesico; *Il Dentista Moderno*, 2013, Novembre: 82-96
- Sisti A., Mottola M.P., Mottola P.;Riabilitazione bilaterale con chirurgia guidata; *Numeri Uno*, 2013; 16:16-18
- Ponzi A.; Echoplan: accuracy dell'implantologia guidata; *Numeri Uno*, 2013;16:12-13
- Morandini E.; La precisione nel cr.co. laser sinterizzato rivestito in ceramica parte 1; *NumeriUno*, 2013; 17: 9-11
- Figliuzzi M. M., De Fazio R., Tiano R., Scordamaglia F., Fortunato L.; Riabilitazione con impianto post-estrattivo immediato in zona estetica: case report; *Numeri Uno*, 17, 2013, 21-22
- Canullo L., Cicchese P., Marinotti F.;Riabilitazione implanto-supportata di entrambi i mascellari edentuli con carico immediato; *Numeri Uno*, 2013; 16, 14-15
- Beolchini M., Lang N.L., Viganò P., Bengazi F., Triana B.G., Botticelli D.; The edentulous ridge expansion (ere) technique an experimental study in the dog; *Clinical Oral Implant Research*, 2013; 25:1207-1211; doi: 10.1111/clr.12263. Epub 2013 Sep 12
- Bengazi F., Botticelli D., Favero V., Perini A., Urbizo Velez J., Lang N.P.; Influence of presence or absence of keratinized mucosa on the alveolar bony crest level as it relates to different buccal marginal bone thicknesses. an experimental study in dogs; *Clinical Oral Implant Research*, 2014; 25:1065-71 ;doi: 10.1111/clr.12233. Epub 2013 Jul 29.
- Crespi R., Capparè P., Gherlone E.F.; Electrical mallet in implants placed in fresh extraction sockets with simultaneous osteotome sinus floor elevation; *The International Journal of Oral & Maxillofacial Implants* 2013;28:869-874; doi: 10.11607/jomi.2679
- Crespi R., Capparè P., Gherlone E.F.; Electrical mallet provides essential advantages in split-crest and immediate implant placement; *International Journal of Oral and Maxillofacial Surgery* 2014;18:59-64; doi: 10.1007/s10006-013-0389-2. Epub 2013 Jan 18
- Csonka M.; Split crest di una cresta molto sottile con il magnetic mallet; *Numeri Uno*, 2013,16:22-23
- Calesini G., Scipioni A.; Approccio rigenerativo sistematico finalizzato all'integrazione morfo-funzionale in implantoprotesi *Numeri Uno*, 16: 6-9, 2013
- Bressan E., Lang N.P., Corazza B., Rizzi S., Almagro Urrutia Z., Botticelli D.; The platform switching concept revisited. an experimental study in dogs. *Clinical Oral Implant Research*, 2013; 00:1-7; doi: 10.1111/clr.12262
- Corrente G., Abundo R., Greppi M., Perelli M., Villa A.; Posizionamento e ricostruzione dei tessuti duri e molli: un protocollo semplificato ; *Numeri Uno*, 2013, 17:14-17
- Guidi R., Viscioni A., Dattola F., Carinci F.; Dental implants inserted in native bone: cases series analyses.; *Dental Research Journal*

- 2012;9:s175-180; doi: 10.4103/1735-3327.109747
- Canullo L., Cicchese P., Marinotti F.; Valutazione di una procedura clinica e tecnica per la riabilitazione di mascellari edentuli; *Il dentista moderno*, 2012; Marzo: 86-102
- Covani U., Ricci M., Tonelli P., Barone A.; An evaluation of new designs in implant-abutment connections: a finite element method assessment; *Implant Dentistry / volume 22, Number 3 2013*; DOI: 10.1097/ID.0b013e318292625f
- Micarelli C., Canullo L., Grusovin M.G., Peñarrocha Oltra D.; Cell adhesion to titanium abutments after different cleaning procedures; *Clinical Oral Implants Research* 2013;24,79-102 (Suppl. 9); doi: 10.11607/jomi.2664
- Canullo L., Peñarrocha Oltra D., Covani U., Micarelli C., Massidda O.; Hard tissue response to plasma of argon cleaning treatment on titanium abutments: 2-year follow-up rct; *Clinical Oral Implants Research* 2013; 24:27-47 (Suppl. 9)
- Canullo L., Peñarrocha Oltra D., Micarelli C., Massidda O., Bazzoli M.; Risposta dei tessuti duri alla pulizia con plasma di argon/sterilizzazione di pilastri in titanio individualizzati, vs pulizia di 5 secondi con vapore: risultati di un studio controllato randomizzato in pazienti con una situazione parodontale favorevole con follow-up a 2 anni dal carico; *European Journal of Oral Implantology* 2013;6(3):251-60
- Canullo L., Peñarrocha Oltra D., Clementini M., Iannello G., Micarelli C.; Impact of plasma of argon cleaning treatment on implant abutments in patients with a history of periodontal disease and thin biotype: radiographic results at 24-month follow-up of a rct; *Clinical Oral Implants Research* 2015;26(1):8-14; doi: 10.1111/clr.12290. Epub 2013 Nov 6
- Canullo L., Cassinelli C., Götz W., Tarnow D.; Plasma of argon accelerates murine fibroblast adhesion in early stages of titanium disk colonization; *The International Journal of Oral & Maxillofacial Implants* 2013;28(4):957-62; doi: 10.11607/jomi.2664
- Avellino W., De Maria A., Milan U., Tamagnone L., Delle Rose D.; Direct prosthetic framework (D.P.F.)
Numeri Uno, 2013; 17:18-20
- Agustín Panadero R., Fons Font A., Román Rodríguez J.L., Solá Ruíz M.F., Cebriá J.R.; Sobredentadura implantosoportada de inserción horizontal; *Gaceta Dental* 249, 2013; 100-112
- Sandri L.P.; Preparazione protesica mediante tecnica B.O.P.T.: caso clinico; *Numeri Uno*, 2013;17:6-8
- Canullo L., Cicchese P., Marinotti F., Sisti A.; Strategia protesica minimamente invasiva negli impianti post-estrattivi: posizionamento e avvistamento; *Il Dentista Moderno*, 2011, Dicembre: 46-54
- Bengazi F., Lang NP, Caroprese M, Velez JU, Favero V, Botticelli D; Dimensional changes in soft tissues around dental implants following free gingival grafting: an experimental study in dogs; *Clinical Oral Implant Research* 26, 176-182, 2015, doi: 10.1111/clr.12280
- Micarelli C., Canullo L., Giuliano I.; Implant/abutment connection deformation after prosthetics procedures - an in vitro study ; *International Journal of Prosthodontics*, 1-9,2014, Early view in ahead of print, accepted July 21st, 2015 doi to be attributed
- Peñarrocha-Oltra D, Covani U, Peñarrocha M, Peñarrocha-Diago M.; Immediate versus conventional loading with fixed full-arch prostheses in mandibles with failing dentition: a prospective controlled study; *International Journal of Oral and Maxillofacial Implants* 30, 2015:427-434; doi: 10.11607/jomi.3534
- Prati C, Zamparini F, Ciulla A, Buonavoglia A, Gatto MR, Piattelli A, Gandolfi MG; Evaluation of marginal bone level of Premium implants; XXIII Congress SIO, Milano 6-7 febbraio Poster; 2015
- Gandolfi MG, Siboni F, Piattelli A, Prati C; Nano-topography, microchemical properties and calcium phosphate nucleation of Premium implants; 30th Annual Congress, American Academy of Osseointegration, San Francisco, 12-14 March Poster , 2015 Id 2088727
- Guazzotti PP; Carico immediato di impianti post estrattivi: presentazione di un caso clinico full-arch; *Doctor Os*, XXVI, 01, gennaio 24-29 ; 2015
- Penarrocha-Oltra D, Rossetti PHO, Covani U, Galluccio F, Canullo L; Microbial leakage at the implant/abutment connection due to implant insertion maneuvers: cross-sectional study 5 years post loading in healthy patients.; *Journal of Oral Implantology*, accepted for publication January 2015
- Agustín-Panadero R., Serra-Pastor B., Chust-López C., Fons-Font A., Ferreiroa A. ; Immediate placement of single implant simultaneously with immediate loading in a fresh socket associated to periapical infection: A clinical case report; *Journal of Clinical and Experimental Dentistry* .;7(1), 2015:175-9
- Canullo L., Peñarrocha-Oltra D., Covani U., Botticelli D., Serino G., Peñarrocha M.; Clinical and microbiological findings in patients with peri-implantitis: a cross-sectional study; *Clinical Oral Implant Research*, 00, 1-7,2015 ; doi: 10.1111/clr.12557
- Mainetti T, Lang NP, Bengazi F, Favero V, Soto Cantero L, Botticelli D; Sequential healing at implants installed immediately into extraction sockets. An experimental study in dogs; *Clinical Oral Implant Research*, 00, 1-9, 2014, doi: 10.1111/clr.12533
- Beolchini M, Lang NP, Gomez Moreno G, Iezzi G, Botticelli D, Calvo Guirado JL; Bone healing at implants with different surface configuration: an experimental study in dogs; *Clinical Oral Implant Research* 00, 1-7, 2015, doi: 10.1111/clr.12562
- Borgia V, Alfonsi F, Toti P, Tonelli P, Covani U, Barone A; Immediate restoration of post-extraction implants. A 7 years prospective single cohort study.; 30th Annual Congress, American Academy of Osseointegration, San Francisco, 12-14 March Poster ; 2015
- Kern JS, Kern T, Wolfart S, Heussen N; A systematic review and meta-analysis of removable and fixed implant-supported prostheses in edentulous jaws: post-loading implant loss; *Clinical Oral Implant Research*, 00, 1-22, 2015, doi: 10.1111/clr.12531
- Crespi R, Bruschi GB, Gastaldi G, Capparè P, Gherlone EF ; Immediate loaded implants in split-crest procedure; *Clin Implant Dent Relat Res.*, Mar 17, 2015 doi: 10.1111/cid.12316
- Martín Anciburo Miguel Ángel; Rehabilitación unitaria implantosoportada utilizando la técnica B.O.P.T.; *Numeri Uno* 04, 2015: 11-14
- Requena Gómez E., Cervantes Haro MN, Aragonese Lamas JM ; ¿Es la cirugía guiada junto a la carga inmediata una técnica predecible? A propósito de un caso clínico; *Numeri Uno* 04, 2015: 16-19
- Canullo L, Peñarrocha-Oltra D, Marchionni S, Bagán L, Peñarrocha-Diago MA, Micarelli C.; Soft tissue cell adhesion to titanium abutments after different cleaning procedures: Preliminary results of a randomized clinical trial.; *Medicina Oral y Patología Oral Cirugía Bucal*, published on line 2013 Oct 13, 2014 Mar 1;19(2): el 77-83, doi: 10.4317/medoral.19329
- Pellicer-Chover H, Peñarrocha-Oltra D, Bagán L, Fichy-Fernandez AJ, Canullo L, Peñarrocha-Diago M; Single blind randomized clinical trial to evaluate clinical and radiological outcomes after one year of immediate versus delayed implant placement supporting full-arch prosthesis; *Medicina Oral y Patología Oral Cirugía Bucal*, 1; 19(3), 2014: 295-301, doi: 10.4317/medoral.19536
- Crespi R, Capparè P, Polizzi E, Gherlone E; Fresh-socket implants of different collar length: Clinical evaluation in the aesthetic zone; *Clinical Implant Dentistry and Related research*, 00, 2014 : 1-8, early view in ahead of print, first published on line 7 Feb 2014 doi 10.1111/cid.12202
- Negri B, López Marí M, Maté Sánchez de Val JE, Iezzi G, Bravo González LA, Calvo Guirado JL; Biological width formation to immediate implants placed at different levels in relation to the crestal bone - an experimental study in dogs; *Clinical Oral Implant Research*, 00, 2014: 1-11, Early view in ahead of print, accepted 06 January 2014 doi 10.1111/clr.12345,
- Strietzel FP, Neumann K, Hertel M ; Impact of platform switching on marginal peri-implant bone-level changes. A systematic review and meta-analysis; *Clinical Oral Implant Research*, 00, 2014: 1-16, Early view in ahead of print, accepted 11 December 2013, doi 10.1111/clr.123339
- Peñarrocha-Oltra D, Candel-Marti E, Peñarrocha-Diago M, Agustín-Panadero R, Canullo L, Peñarrocha M; The Horizontal Denture©: a prosthodontic alternative for Severe Maxillary Atrophy. A technical note; *Journal of Oral Implantology*, Early view in ahead of print, accepted 8 January 2014, 2014
- Maiorana C, Farronato D, Pieroni S, Cicciù M, Andreoni D, Santoro F; A four-year survival rate multicenter prospective clinical study on 377 implants - correlations between implant insertion torque, diameter and bone quality; *Journal of Oral Implantology*, 2014, Early view in ahead of print, accepted 11 February 2014
- Canullo L, Peñarrocha-Oltra D, Soldini C, Mazzocco F, Peñarrocha M, Covani U; Microbiological assessment of the implant-abutment interface in different connections: cross-sectional study after 5 years of functional loading; *Clinical Oral Implant Research*, 00, 2014: 1-9, Early view in ahead of print, accepted 22 February 2014, doi 10.1111/clr.12383
- Mainetti T, Lang N, Bengazi F, Sbricoli L, Soto Cantero L, Botticelli D.; Immediate loading of implants installed in a healed alveolar bony ridge or immediately after tooth extraction: an experimental study in dogs; *Clinical Oral Implant Research*, 00, 2014: 1-8, Early view in ahead of print, accepted 5 March 2014, doi 10.1111/clr.12389
- Engelhardt S, Papacosta S, Rathe F, Ozen J, Jansen J.A., Junker R.; Annual failure rates and marginal bone-level changes of immediate compared to conventional loading of dental implants. A systematic review of the literature and meta-analysis; *Clinical Oral Implant Research*, 00, 2014: 1-17, Early view in ahead of print, accepted 9 February 2014, doi 10.1111/clr.12363
- Bruschi GB, Crespi R, Capparè P, Grande N, Bruschi E, Gherlone E; Radiographic evaluation of crestal bone levels of delayed implants at 5 medium term follow up; *International Journal of Oral & Maxillofacial Implants*, 29;2014: 441-447 doi 10.11607/jomi.3254
- Sbordone C, Toti P, Martuscelli R, Guidetti F, Sbordone L, Ramaglia L; A 5-year implant follow-up in maxillary and mandibular horizontal osseous onlay grafts and native bone; *Journal of Oral Implantology*,

- Early view in ahead of print, accepted 4 March 2014; 2014
- Canullo L, Micarelli C, Bettazzoni L, Magnelli A, Baldissara P; Shear bond strength of veneering porcelain to zirconia after argon plasma treatment; *International Journal of Prosthodontics*, Mar-Apr, 27(2), 2014: 137-9, 2014 doi: 10.11607/ijp.3722
- Canullo L, Micarelli C, Bettazzoni L, Koçi B, Baldissara P; Zirconia-Composite bonding after plasma of argon treatment; *International Journal of Prosthodontics*, 27:267-269, 2014, doi: 10.11607/ijp.3686
- Peñarrocha-Oltra D, Peñarrocha-Diago M, Canullo L, Covani U, Peñarrocha Miguel; Patient-reported outcomes of immediate versus conventional loading with fixed full-arch prostheses in the maxilla: a non-randomized controlled prospective study; *The International Journal of Oral & Maxillofacial Implants*, 29 (3), 690-698; 2014
- Covani U, Canullo L, Toti P, Alfonsi F, Barone A; Tissue stability of implants placed in fresh extraction sockets - a 5 year prospective single cohort study; *Journal of Periodontology*, 85: 323-332, 2014, doi 10.1902/jop2014.140175
- De Santis E, Lang NP, Favero G, Beolchini M, Morelli F, Botticelli D.; Healing at mandibular block-grafted sites. An experimental study in dogs; *Clinical Oral Implant Research*, 00, 2014: 1-7, Early view in ahead of print, accepted 17 May 2014, doi 10.1111/clr.12434
- Crespi R, Bruschi GB, Capparè P, Gherlone E.; The utility of the electric mallet; *The Journal of Craniofacial Surgery*, 25 May (3), 793-795, 2014, doi 10.1097/SCS.0000000000000523;2014
- Covani U, Marconcini S, Ferrini F, Gelpi F, Finotti M, Barone A.; Post-traumatic use of dental implants immediately after tooth extraction - clinical study; *The Journal of Craniofacial Surgery*, 25 May (3), 796-798, 2014, doi 10.1097/SCS.0000000000000522
- Calesini G, Zarone F, Sorrentino R, Micarelli C, Fabianelli A, Paccchini F, Gherlone E.; Effect of 2 impression techniques on the dimensional accuracy of working implant prosthesis models - an in vitro study; *The Journal of Craniofacial Surgery*, 25 May (3), 822-827, 2014, doi 10.1097/SCS.0000000000000715
- Quaranta A, Andreana S, Pompa G, Procaccini M; Active implant peri-apical lesion - a case report treated via guided bone regeneration with a 5-year clinical and radiographic follow-up; *Journal of Oral Implantology*, 40 (3), 313-319, 2014, doi: 10.1563/AAI.JOI.D.11.00214
- J. Viña-Almunia ; Microbial colonization of the implant connection with cemented versus screw-retained suprastructures; Oral presentation ,Clinical research - Prosthodontically oriented
- EAO Congress, Rome 25-27 September 2014, *Clinical Oral Implant Research*, 25 (suppl. 10), 93, 2014
- Cocchetto R.; Improved cementation technique for implant restorations to avoid periimplant cement remnants: clinical and microscopical evaluation with two different abutment design; Oral presentation, Clinical research - Prosthodontically oriented
- EAO Congress, Rome 25-27 September 2014, *Clinical Oral Implant Research*, 25 (suppl. 10), 96; 2014
- Augustín-Panadero R, Serra-Pastor B, Roig-Vanaclocha A, Román-Rodríguez JL, Fons-Font A; Mechanical behavior of provisional implant prosthetic abutments; *Medicina Oral Patología Oral y Cirugía Bucal*, 1-9, 2014, Early view in ahead of print, accepted July 2014, doi 10.4317/medoral.19958.
- Micarelli C, Canullo L, Giuliano I.; Implant/abutment connection deformation after prosthodontics procedures - an in vitro study ; *International Journal of Prosthodontics*, 1-9,2014, Early view in ahead of print, accepted July 21st, 2015 doi to be attributed
- Kinaia BM, Shah M, Neely AL, Goodies HE; Crestal bone level changes around immediately placed implants - A systematic review and meta-analyses with at least 12 months follow up after functional loading; *Journal of Periodontology*, 2014, early view in ahead of print, doi: 10.1902/jop2014.130722;2014
- Cicciù M, Bramanti E, Matacena G, Guglielmino E, Risistano G.; FEM evaluation of cemented-retained versus screw-retained dental implant single-tooth crown prosthesis; *International Journal of Clinical and Experimental Medicine* 7(4) , 2014: 817-825; doi: 1940-5901. ijcem.1402025
- Crespi R, Capparè P, Polizzi EM, Gherlone EF; Tissue remodeling after bone expansion in grafted and ungrafted sockets; *The International Journal of Oral and Maxillofacial Implants*, 29, 2014: 699-704, , doi: 10.11607/jomi.3535
- Bruschi GB, Crespi R, Capparè P, Gherlone E.; Clinical Study of flap design to increase the keratinized gingiva around implants - a 4 year follow-up; *Journal of Oral Implantology*, 40(4) , 2014: 459-464, doi: 10.1563/aaid-joi-d-11-00236
- Romanos GR, Javed F.; Platform switching minimises crestal bone loss around dental implants - truth or myth?; *Journal of Oral Rehabilitation*, 2014, early view in ahead of printing, accepted for publication 30 April 2014, doi 10.1111/joor.12189
- Gaspari L.; Tecnica conometrica con provvisorio elettrosaldato per carico immediato; *Italian Dental Journal*, agosto, 29; 2014
- Lumetti S, Galli C, Smerieri A, Macaluso G, Manfredi E, Ghiacci G, Di Blasio A, Megri M.; The effect of age, gender and insertion site on marginal bone loss around endosseous implants: results for a 3 year trial; Poster, EAO Congress, Rome 25-27 September 2014, *Clinical Oral Implant Research*, 25 (suppl. 10), 440; 2014
- Lumetti S, Di Blasio A, Manfredi E, Ghiacci G, Toffoli A, Bonanini M, Macaluso G, Galli C.; Implant surface microtopography affects the pater of cell growth, cell-to-cell contacts and the expression of Connexin 43;Poster, EAO Congress, Rome 25-27 September 2014, *Clinical Oral Implant Research*, 25 (suppl. 10), 222; 2014
- Caneva M, Lang NP, Calvo Guirado JL, Spriano AM, Iezzi G, Botticelli D.; Bone healing at bicortically installed implants with different surface configurations. An experimental study in rabbits; *Clinical Oral Implant Research*, 00, 2014: 1-7, Early view in ahead of printing, accepted 29 July 2014, doi:10.1111/clr.12475
- D'Ercole S, Tripodi D, Marzo G, Bernardi S, Continenza MA, Piatelli A, Iaculli F, Mummolo S.; Microleakage of bacteria in different implant-abutment assemblies: an in vitro study; *Journal of Applied Biomaterial and Functional Materials*, 2014, accepted June 12, 2014, doi: 105301/jabfm.5000214
- Peñarrocha-Oltra D, Peñarrocha-Diago M, Aloy-Prosper A, Covani U, Peñarrocha M.; Immediate versus conventional loading of complete-arch implant-supported prostheses in mandibles with failing dentition: a patient centered controlled prospective study; *Journal of oral and Maxillofacial Implants*, submitted; 2014
- Bowen Antolin A, Ariño B, Arlandi Garrido M.; Regeneración ósea periimplantaria con fosfato de calcio bifásico y ácido poliláctico; *Gaceta Dental*, 260(7), 2014: 174-186;
- Morelli F, Lang NP, Bengazi F, Baffone D, Vila Morales CD, Botticelli D.; Influence of bone marrow on osseointegration in long bones: an experimental study in sheep; *Clinical Oral Implant Research*, 00, 1-7, 2014, Early view in ahead of printing, accepted 29 August 2014, doi:10.1111/clr.12487
- Marchetti E, Ratta S, Mummolo S, Tecco S, Pecci R, Bedini R, Marzo G.; Evaluation of an endosseous oral implant system according to UNI EN ISO 14801 Fatigue Test Protocol; *Implant Dentistry*, 2014, Early View in ahead of print, doi: 10.1097/id.151
- Negri M, Galli C, Smerieri A, Macaluso GM, Manfredi E, Ghiacci G, Toffoli A, Bonanini M, Lumetti S; The effect of age, gender and insertion site on marginal bone loss around endosseous implants: results from a 3-year trial with Premium Implant System; *BioMed research International*, Volume 2014, Article ID 369051, 7 pages, doi. org/10.1155/2014/369051
- Esposito M, Ardebili Y, Worthington HV; Interventions for replacing missing teeth: different types of dental implants (Review); *The Cochrane Collaboration*, John Wiley and Sons, Ltd; 2014
- Mainetti T, Lang NP, Bengazi F, Favero V, Soto Cantero L, Botticelli D; Sequential healing at implants installed immediately into extraction sockets. An experimental study in dogs; *Clinical Oral Implant Research*, 00, 1-9, 2014, doi: 10.1111/clr.12533
- Crespi R, Capparè P, Gastaldi G, Gherlone EF; Immediate Occlusal loading of full-arch rehabilitations: screw-retained versus cement-retained prosthesis. An 8 year clinical evaluation; *International Journal of Oral & Maxillofacial Implants* 29, 2014:1406-1411; doi: 10.11607/jomi.3746
- Pradies Ramiro G., Abad Coronel C., García Martínez I., Ferreiroa Navarro A.; Impresiones fiables: dos propuestas para un mismo objetivo; *Numeri Uno*, 01, 2014, 6-9
- Machín Muñoz A.; Regeneración ósea y gingival en implantes inmediatos post-extracción; *Numeri Uno* 01 , 2014: 20-21
- Loi I.; Riabilitazione implanto-protetica di elemento incisivo frontale con impianto Prama; *Numeri Uno* 20, 2014: 12-13
- Loi I.; Riabilitazione implanto-protetica di ponte distale con impianti Prama; *Numeri Uno* 20, 2014: 14-15
- Gorni F.; Riabilitazione di elemento singolo in zona estetica con impianto Prama RF; *Numeri Uno* 20, 2014: 16-17
- Vedove F.; Riabilitazione di elemento singolo in zona estetica con impianto Prama; *Numeri Uno* 20, 2014: 18-19
- Andreoni D.; Riabilitazione di elemento singolo in posizione 4.6 con impianto Prama; *Numeri Uno* 20, 2014: 20-21
- Sandri L.P.; Utilizzo clinico dei nuovi impianti Prama: inserimento e riabilitazione con un singolo impianto; *Numeri Uno* 20, 2014: 22-24
- Csonka M.; Trattamento implantologico delle creste sottili: Split Crest o GBR?; *Numeri Uno* 19, 2014: 12-14
- Vischia F., Roncoroni F.; Ortodonzia protesica mediante tecnica B.O.P.T.; *Numeri Uno* 19, 2014: 19-21
- Schirripa G., Schirripa F.; Carico immediato; *Numeri Uno* 19, 2014: 22-24
- Baldi D., Colombo J., Pera P., Hauschild U.; Una tecnica minimamente invasiva: implantologia con utilizzo di impianti a diametro ridotto e tecniche CAD CAM per una provvisorizzazione a lungo termine; *Numeri Uno* 18, 2014: 6-9
- Gaspari L.; Implantoprotesi conometrica elettrosaldato chairside a carico immediato - caso clinico; *Numeri Uno*, 18, 2014:12-14
- Loi I.; Tecnica B.O.P.T. su denti e impianti per la riabilitazione di un'ar-

- cata completa; Numeri Uno 18, 2014:21-22
- Morandini E.; La precisione nel Cr.Co. laser sinterizzato rivestito in ceramica; NumeriUno 17, 2013: 9-11 - NumeriUno 18, 2014: 16-19
- Loi I.; Técnica B.O.P.T. sobre dientes naturales; Numeri Uno 02, 2014: 8-9
- Loi I.; Técnica B.O.P.T. sobre dientes e implantes para la rehabilitación de los dos arcos completos; Numeri Uno 02, 2014 : 14
- Canullo L, Cassinelli C, Goetz W, Tarnow D; Il plasma di argon accelera l'adesione dei fibroblasti murini nelle fasi precoci della colonizzazione di dischetti in titanio; International Journal of Oral and Maxillofacial Implants 2013; 28: 957-962. DOI: 10.11607/jomi.2664
- Bengazi F, Botticelli D, Favero V, Perini A, Urbizo Velez J, Lang NP; Influence of presence or absence of keratinized mucosa on the alveolar bony crest level as it relates to different buccal margin bone thicknesses. An experimental study in dogs; Clinical Oral Implant Research, 00, 2013, 1-7, Accepted 26 June 2013, first published online on 29/07/2013, DOI 10.1111/clr.12233
- Peñarrocha-Oltra D, Covani U, Aparicio A, Ata-Ali J, Peñarrocha-Diogo Miguel, Peñarrocha-Diogo María; Immediate versus conventional loading for the maxilla with implants placed into fresh and healed extraction sites to support a full-arch fixed prosthesis: nonrandomized controlled clinical study; International Journal of Oral and Maxillofacial Implants 2013; 28: 1116-1124 DOI: 10.11607/jomi.3119
- Covani U, Ricci M, Tonelli P, Barone A; An evaluation of new designs in implant-abutment connections: a finite element method assessment; Implant Dentistry, 2013, Jun22(3): 263-267, DOI 10.1097/ID.0b013e318292625f
- Crespi R, Capparè P, Gherlone EF; Electrical mallet in implants placed in fresh extraction sockets with simultaneous osteotome sinus floor elevation; International Journal of Oral and Maxillofacial Implants, 2013; 28(3): 869-874, doi: 10.11607/jomi.2679
- Panadero RA, Fons Font A, Granell Ruiz M, Román Rodríguez JL, Solá Ruiz MF, Rubio Cebriá J; Sobredentadura implantosoportada de inserción horizontal; Gaceta Dental, 249: 100-112, 2013
- Beolchini M, Lang NL, Viganò P, Bengazi F, Triana BG, Botticelli D; The edentulous ridge expansion (ERE) technique - an experimental study in dogs; Clinical Oral Implant research, 2013: 1-7, published online early view in ahead of print in September 2013, doi: 10.1111/clr.12262
- Bressan E., Lang NP, Corazza B, Rizzi S, Almagro Urrutia Z, Botticelli D; The Platform Switching concept revisited. An experimental study in dogs; Clinical Oral Implant research, 2013: 1-5, published online early view in ahead of print in September 2013, doi: 10.1111/clr.12263
- Crespi R, Capparè P, Gherlone EF; Electrical mallet provides essential advantages in split-crest and immediate implant placement; Oral and Maxillofacial Surgery, 2013, (18): published online early view in ahead of print in January 2013, doi: 10.1007/s10006-013-0389-2
- Canullo L, Peñarrocha-Oltra D, Marchionni S, Bagán L, Peñarrocha-Diogo MA, Micarelli C.; Soft tissue cell adhesion to titanium abutments after different cleaning procedures: Preliminary results of a randomized clinical trial.; Medicina Oral y Patología Oral Cirugía Bucal, published online 2013 Oct 13, 2014 Mar 1;19(2): el 77-83, doi: 10.4317/medoral.19329
- Canullo L, Peñarrocha D, Peñarrocha M, Rocío A-G, Peñarrocha-Diogo M.; Piezoelectric vs. conventional drilling in implant site preparation: pilot controlled randomized clinical trial with crossover design.; Clinical Oral Implant Research 00, 2013, 1-8, published online early view in ahead of print in October 2013, doi: 10.1111/clr.12278
- Micarelli C, Canullo L, Grusovin MG, Peñarrocha Oltra D.; Cell adhesion to titanium abutments after different cleaning procedures; Clinical Oral Implant Research, 24(Suppl.9), 2013 : 79-102
- Canullo L, Peñarrocha D, Covani U, Micarelli C, Massidda O.; Hard Tissue response to plasma of argon cleaning treatment on titanium abutments - 2 year follow-up RCT; Clinical Oral Implant Research, 24(Suppl.9), 27-47, 2013
- De Risi V, Clementini M, Vittorini G, Mannocci A, De Sanctis M; Alveolar ridge preservation techniques: a systematic review and meta-analysis of histological and histomorphometrical data; Clinical Oral Implant Research, 00, 2013: 000-000, Early view in ahead of print, accepted September 2013, doi 10.1111/clr.12288
- Canullo L, Peñarrocha D, Clementini M, Iannello G, Micarelli C; Impact of plasma of argon cleaning treatment on implant abutments in patients with a history of periodontal disease and thin biotype - radiographic results at 24 months follow-up of a RCT; Clinical Oral Implant Research, 00, 2013: 000-000, Early view in ahead of print, accepted 18 September 2013, doi 10.1111/clr.12290
- Canullo L, Peñarrocha D, Micarelli C, Massidda O, Bazzoli M; Hard tissue response to argon plasma cleaning / sterilization of customized titanium abutments versus 5-second steam cleaning: results of a 2-year post-loading follow-up from an explanatory randomized controlled trial in periodontally healthy patients; European Journal of Oral Implantology, Autumn ; 6(3) ,2013:251-60
- Petrillo N.; Carico immediato full arch mascellare e mandibolare: un nuovo approccio chirurgico e protesico; Il Dentista Moderno, 2013 Novembre 2013: 82-96
- Baffone G, Lang NP, Pantani F, Favero G, Ferri M, Botticelli D; Hard and soft tissue changes around implants installed in regular-sized and reduced alveolar bony ridges. An experimental study in dogs; Clinical Oral Implant Research, 00, Early view in ahead of print, accepted 28 October 2013: 1-6, doi 10.1111/clr.12306
- Beolchini M, Lang NL, Ricci E, Bengazi F, Garcia Triana B, Botticelli D; Influence on alveolar resorption of the buccal bony plate width in the edentulous ridge expansion (E.R.E.) - an experimental study in the dog; Clinical Oral Implant Research, 00, 2013: 1-6, Early view in ahead of print, accepted 28 October 2013doi 10.1111/clr.12308
- Strietzel FP, Neumann K, Hertel M; Impact of platform switching on marginal peri-implant bone-level changes. A systematic review and meta-analysis; Clinical Oral Implant Research, 00, 2014: 1-16, Early view in ahead of print, accepted 11 December 2013, doi 10.1111/clr.123339
- Morandini E.; La precisione nel Cr.Co. laser sinterizzato rivestito in ceramica; NumeriUno 17, 2013: 9-11 - NumeriUno 18, 2014: 16-19
- Sandri L.P.; Preparazione protesica mediante tecnica B.O.P.T.: caso clinico; Numeri Uno 17, 2013 :6-8
- Corrente G., Abundo R., Greppi M., Perelli M., Villa A.; Posizionamento e ricostruzione dei tessuti duri e molli: un protocollo semplificato; Numeri Uno 17, 2013:14-17
- Avellino W., De Maria A., Milan U., Tamagnone L., Delle Rose D.; Direct Prosthetic Framework (D.P.F.); Numeri Uno, 17, 2013: 18-20
- Figliuzzi M. M., De Fazio R., Tiano R., Scordamaglia F., Fortunato L.; Riabilitazione con impianto post-estrattivo immediato in zona estetica: Case Report; Numeri Uno 17, 2013:21-22
- Fadda M.; Caso clinico con M.F. Extrusion; Numeri Uno, 17, 2013:26
- Cardarelli F.; Effetti dentofacciali della terapia ortodontica in dentizione mista per la correzione delle II Classi; Numeri Uno 17, 2013: 28-31
- Calesini G., Scipioni A.; Approccio rigenerativo sistematico finalizzato all'integrazione morfo-funzionale in implantoprotesi; Numeri Uno 16, 2013: 6-9
- Ponzi A.; Echo Plan: accuracy dell'implantologia guidata; Numeri Uno 16, 2013: 12-13
- Canullo L., Cicchese P., Marinotti F.; Riabilitazione impianto-supportata di entrambi i mascellari edentuli con carico immediato; Numeri Uno 16, 2013: 14-15
- Sisti A., Mottola M.P., Mottola P.; Riabilitazione bilaterale con chirurgia guidata; Numeri Uno 16, 2013: 16-18
- Csonka M.; Split crest di una cresta molto sottile con il Magnetic Mallet; Numeri Uno 16, 2013: 22-23
- Guidi R, Viscioni A, Dattola F, Carinci F; Dental implants inserted in native bone: cases series analyses; Dental Research Journal, 12(9), Issue 8 (Suppl Issue 2), 175-180; 2012
- Canullo L, Cicchese P, Marinotti F.; Valutazione di una procedura clinica e tecnica per la riabilitazione dei mascellari edentuli; Il Dentista Moderno, Marzo: 86-102, 2012
- Canullo L, Cicchese P, Marinotti F, Sisti A; Strategia protesica minimamente invasiva negli impianti post-estrattivi: posizionamento e avvitamento immediato del moncone definitivo per protesi su impianti secondo il concetto del platform switching; Il Dentista Moderno, Dicembre: 46-54, 2011
- Segura Andrés G., Martínez Lage J. F., Ferreira A., Faus López J., Agustín Panadero R.; Rehabilitación protesica en un maxilar atrófico a consecuencia de un trauma facial; Gaceta Dental, 2013; 244:112-118
- Canullo L., Micarelli C., Clementini M.; Hard tissue response to argon plasma cleaning treatment on titanium abutments: 2-year follow-up RCT; European Journal of Oral Implantology, 6 (Suppl. Spring), 2013: S21-S22
- Rossi F., Lang N. P., De Santis E., Morelli F., Favero G., Botticelli D.; Bone-healing pattern at the surface of titanium implants: an experimental study in the dog; Clinical Oral Implant Research, Early View, First Published online on 2013, January 4th; DOI: 10.1902/jop.2010.100428
- Clementini M., Canullo L., Micarelli C.; Fibroblast growth on titanium disks treated by argon plasma: an in vitro triple-blinded study; European Journal of Oral Implantology, 6 (Suppl. Spring), 2013: S29-S30
- Canullo L., Micarelli C., Iannello G.; Microscopical and chemical surface characterization of the gingival portion and connection of an internal hexagon abutment before and after different technical stages of preparation; Clinical Oral Implant Research, 2013, 24: 606-611; DOI: 10.1111/j.1600-0501.2012.02499.x
- Canullo L., Heinemann F., Gedrange T., Biffar R., Kunert-Keil C.; Histological evaluation at different times after augmentation of extraction sites grafted with a magnesium-enriched hydroxyapatite: double-blinded randomized controlled trial; Clinical Oral Implant Research, Early View, First Published online on 2013, January 4th; DOI: 10.1111/clr.12035



rev.07-18



Sweden & Martina S.p.A.

Via Veneto, 10
35020 Due Carrare (PD), Italy
Tel. +39.049.9124300
Fax +39.049.9124290
info@sweden-martina.com

www.sweden-martina.com

Sweden & Martina Ltd

Unit 45
Basepoint Business Centre
Metcalf Way, Crawley,
West Sussex, RH11 7XX, UK
Toll free 0800 1123575
info.uk@sweden-martina.com

Sweden & Martina Inc.

Distributor for U.S.
80 Broad Street, 5th floor, Suite 502
New York, NY 10004
Toll free 1-844-8MARTINA(1-844-862-7846)
info.us@sweden-martina.com

www.sweden-martinainc.com

Sweden & Martina Mediterranea S.L. - España

info.es@sweden-martina.com
Sweden & Martina Lda - Portugal
info.pt@sweden-martina.com

The implants, standard prosthetic components and surgical instruments contained in this catalogue are Medical devices and are manufactured by Sweden & Martina S.p.A. They conform to the ISO 9001 and ISO 13485 standards and are certified with the CE Mark (Class I) and CE 0476 mark (Class IIA and class IIB) in compliance with Regulation (EU) Medical Devices n. 2017/745.

We have met the good manufacturing standards (GMP) set forth by many countries worldwide, including the United States FDA.



Some products may not be regulatory/released for sale in all markets.

All trademarks herein are the property of Sweden & Martina S.p.A. unless otherwise indicated.

This material is intended for laboratories and clinicians and is not intended for patient distribution.

This material is not to be redistributed, duplicated, or disclosed without the express written consent of Sweden & Martina S.p.A.

For additional product information, including indications, contraindications, warnings, precautions, and potential adverse effects, see Sweden & Martina S.p.A. website.

The contents are updated at the time of publication. Check with the company for any subsequent updates.